







Support for your
decontamination process.








PRODUCT CATALOGUE

CLINIPAK
trusted solutions at every step

Table of contents

	Introduction	03-04		Sterilisation	39-42
	Welcome to Clinipak	03		Cleaning and validation	40
	Global presence	04		Sterilisation containers	41
	Industry Partnerships	04			
	Accreditations and memberships	04		Transport	43-48
	Cleaning	05-08		Transportation boxes and crates	44
	Cleaning brushes for cannulae	06		Transportation accessories	45
	Brushes for wash room cleaning	08		Transportation pouches	48
	Disinfection	09-16		Storage	49-52
	Trays and baskets	10		Storage systems	50
				Transportation trays and bags	51
	Inspection	17-22		Use	53-56
	Absorbent towels	18		Theatre sundries	54
	Identification	18			
	Tracking and tracing	19		Index (A-Z)	59
	Packing	23-38			
	Barrier wraps	24			
	Wrap bags	26			
	Transport wraps	27			
	Autoclave tapes	28			
	Protective packaging	29			
	Absorbent liners	32			
	Bags, pouches & reels	34			

Key to icons

	Suitable for Steam Sterilisation		Single Use
	Suitable for Sterilisation with EO Gas		CE Marked to MDD 93/42
	Suitable for Hydrogen Peroxide Sterilisation		Latex Free
	Suitable for Gamma Sterilisation		For UK Distribution Only
	Suitable for Wash/Decontamination		



Welcome to Clinipak

Clinipak is a medical manufacturing and distribution company that provides a distinct range of solutions to support all stages of the decontamination process.

We are one of the only UK-based companies to specialise exclusively in the provision of Central Sterile Supplies and Theatre Sterile Supplies.

For over 30 years, Clinipak has not only been dedicated to supporting the cleaning, disinfection and sterilisation of reusable invasive medical devices, but has maintained a company ethos which focuses on quality, innovation and customer satisfaction.

Our quality management system is not only accredited under the ISO 9001:2008 certification, but enables us to implement more efficient ways of working, improve cost control and develop new techniques to advance current working practices.


As a result of effectively interacting with our customers and understanding their specific needs, we are constantly improving our existing product portfolio and developing innovative, trusted solutions for all central sterile service departments. Our bespoke solutions have been notably designed to overcome common issues such as wet packs, torn sets and handling issues.

Our location in the South East of England and proximity to London provide us with excellent transport links which ensure a high-quality delivery service throughout the UK.

As a unique proposition to our customers, we can provide training, support and clinical guidance to ensure the most cost effective and reliable solutions for individual needs.

Global presence

Clinipak has developed strong relationships on an international scale over the last 5 years and now actively distributes its products to more than 40 countries, particularly across Europe and North America. Please contact our team if you would be interested in becoming a Clinipak distributor.



Industry partnerships

Since its inception in 1983, Clinipak has proudly started to distribute a wide range of worldwide, reputable brands to more effectively exceed the needs of our customers in the UK and overseas. The products of these key suppliers not only complement our product ranges, but support our intentions to provide the most trusted solutions at every step of the decontamination process without fail.

Clinipak is proudly distributing the following brands in the UK:



Technical Tapes



MEDICA
EUROPE




SAFE & SMART HEALTHCARE SOLUTIONS



RØNVIG
Dental Mfg. A/S

Accreditations and memberships





Cleaning

pages 05-08

Cleaning brushes for cannulae

A range of disposable brushes which have been designed for the manual cleaning of dirty cannulated instruments.

There are various styles and sizes available to provide the most effective cleaning solution.

Flexistem



A unique disposable cleaner for all small cannulae, Clinipak Flexistem has been produced in a wide range of materials and various diameters to provide either soft or hard cleaning properties.



CODE	DESCRIPTION	MATERIAL	QUANTITY
FLEX001	3mm soft cleaner	Cotton	1 x 10m coil
FLEX002	3mm super bristle	Cotton & Nylon	1 x 10m coil
FLEX003	3mm double clean 15cm	Cotton & Nylon	100
FLEX004	9mm fleece stem 30cm	Polyester	50
FLEX005	12mm fleece stem 30cm	Polyester	50
FLEX006	15mm bristle stem per 30cm	Polyester & Nylon	20
FLEX007	15mm bristle stem	Polyester & Nylon	1 x 5m coil
FLEX033	3mm fleece stem 30cm	Polyethylene	50

Wavy flexistem



Wavy flexistem has been developed for cleaning the inside of long instruments and prevents the brush from buckling when a harsher pushing force is required.



CODE	DESCRIPTION	MATERIAL	QUANTITY
FLEX007A	15mm wavy bristle stem	Polyester & Nylon	1 x 5m coil
FLEX007B	15mm wavy bristle stem	Polyester & Nylon	10 x 66.5cm coil

Microclene




Microclene is a solution to facilitate the cleaning of extra small bore items. A development of flexistem, microclene is finer, more brush-like, and can be used to safely dislodge deposits during the washing process.



CODE	DESCRIPTION	MATERIAL	QUANTITY
FLEX021	3mm microclene (super fine)	Polyester & Nylon	1 x 7.5m coil
FLEX022	6mm microclene (medium fine)	Polyester & Nylon	1 x 7.5m coil
FLEX023	3mm microclene (super brush)	Polyester & Nylon	1 x 7.5m coil
FLEX024	6mm microclene (medium brush)	Polyester & Nylon	1 x 7.5m coil

Flexistem and microclene dispensers

Clinipak dispensers ensure that flexistem and microclene are accessible at all times. A cutter is also supplied alongside these dispensers to enable users to cut these brushes to any required length with ease.



FLEXISTEM DISPENSER

CODE	DESCRIPTION
FLEX008	Dispenser for the flexistem range, which includes one cutter and one of each FLEX001 – FLEX007. Size: 285 x 190 x 180mm.
FLEX009	Complete refill pack for flexistem dispenser which contains one of each FLEX001 – FLEX007.



MICROCLENE DISPENSER

CODE	DESCRIPTION
FLEX025	Dispenser for the microclene range which includes one cutter and one of each FLEX021 – FLEX024. Size: 285 x 190 x 180mm.
FLEX026	Complete refill pack for microclene dispenser which contains one of each FLEX021 – FLEX024.

"... We have consistently relied on Clinipak over the last decade for its variety of high-quality products that has helped us to achieve excellence in the delivery of our hospital sterile services."

Val O'Brien
Global Technical Lead HS
Synergy Health Plc



Brushes for wash room cleaning

A range of disposable and reusable cleaning brushes for manual cleaning in the wash room.

Blue brushes



A range of brushes with strong nylon bristles mounted into durable handles which are long lasting and safe to use. These brushes are suitable for cleaning and disinfection, but not for steam sterilisation.



CODE	DESCRIPTION
BBR001	Medium blue cleaning brush
BBR002	Round blue cleaning brush
BBR003	Small blue cleaning brush

Endoscopy brushes



Single-use brush for the safe cleaning of endoscopes with a tapered brush head for optimum cleaning capability.



CODE	DESCRIPTION	QUANTITY
BBR005	Brush head 7.5mm (max diameter) x 9mm.	20 x 2.3m coil

Diathermy cleaning brushes



A robust pair of stainless steel brushes for cleaning the tips of diathermy forceps. The removable guard prevents damage to the insulation of instruments and the brush can be separated for safer cleaning.



CODE	DESCRIPTION
DIA001	Supplied as a pair – 1 x 10mm wide and 1 x 15mm wide



Disinfection

pages 09-16

Trays and baskets

A range of high-quality stainless steel washing and instrument baskets produced in line with all industry standard sizes or made to customer requirements.

Clinipak premium baskets

S EO HP W CE LF Custom sizes available

Clinipak premium baskets are produced with high-quality, perforated stainless steel which guarantees both strength and durability at an attractive price. Each basket is resistant to accidental damage, supplied with or without handles, and designed with rounded corners to avoid torn wraps from sharp edges.



CODE	DESCRIPTION	SIZE
BSK316	Din	480 x 250 x 50mm
BSK317	Half Din	250 x 250 x 50mm
BSK318	Double Din	480 x 480 x 50mm
LID316	Lid to fit BSK316	480 x 250mm
LID317	Lid to fit BSK317	250 x 250mm
LID318	Lid to fit BSK318	480 x 480mm

Flat base economy baskets

S EO HP W CE LF

A selection of two baskets with perforated stainless steel sides and a wire mesh base which facilitates effective drainage in the washer disinfector. All baskets are electrolytically polished and are supplied with inner mounted drop handles. Lids are also available to fit these baskets.



CODE	DESCRIPTION	SIZE
CIL-0234	Flat base basket 240	240 x 250 x 60mm
CIL-0237	Flat base basket 480	480 x 250 x 60mm

Wire mesh economy baskets

S EO HP W CE LF

A selection of two baskets made entirely from a 4.2mm woven stainless steel mesh which exhibits strong water drainage properties. The baskets are electrolytically polished and have a stackable rim with a burr-free smooth finish.



CODE	DESCRIPTION	SIZE
CIL0202	Wire mesh basket 240	240 x 250 x 60mm
CIL0210	Wire mesh basket 480	480 x 250 x 60mm

Tray dividers

S EO HP W LF

Clinipak tray dividers are available as a one-piece six compartment divider or adjustable slotted divider bars which can enable users to adjust layouts. A DIN tray would require one DIV480 and two DIV250 to produce 6 compartments. All dividers are made from Stainless steel.




CODE	DESCRIPTION	SIZE
DIV116	A fitted divider with six compartments.	480 x 250 x 38mm
DIV480	Adjustable slotted divider bar to fit trays 480mm long/wide.	480 x 38mm
DIV250	Adjustable slotted divider bar to fit trays 250mm long/wide.	250 x 38mm

Trinket baskets

S EO HP W LF

Small stainless steel baskets to hold fine items such as screws and bolts.



CODE	DESCRIPTION	SIZE
CIL0406	Fine mesh basket	80 x 80 x 34mm
CIL0705	Trinket box	90 x 60 x 25mm
CIL0706	Trinket box	105 x 70 x 25mm
CIL0707	Trinket box	120 x 80 x 25mm
CIL76-0204	Fine mesh basket	90 x 60 x 25mm
CIL76-0185	Trinket box (wire mesh lid)	105 x 70 x 25mm

Nic-nac baskets

S EO HP W CE LF Custom sizes available

A basket selection designed to transport small items through the wash and autoclave process. All nic-nac baskets are produced with a perforated stainless steel sheet with a scalloped hinged lid.



CODE	SIZE
NIC001R	80 dia x 30mm
NIC002	178 x 127 x 25mm
NIC003	178 x 229 x 25mm
NIC003DEEP	180 x 230 x 50mm
NIC004DEEP	180 x 180 x 50mm
NIC005DEEP	105 x 70 x 25mm

Belintra premium baskets



A high quality range of stainless steel trays and baskets with perforated sides and a mesh base to minimise wet packs. Available in three main industry standard sizes.

High-quality Belintra silicone matting can be supplied for each premium basket size.





CODE	SIZE
BELBSK116	480 x 250 x 60mm
BELBSK117	240 x 250 x 60mm
BELBSK110	120 x 250 x 60mm

Belintra Steri-ID range



The Belintra steri-ID range represents a high-quality alternative to the traditional process for tracking the contents of instrument baskets with barcode stickers due to the use of embedded tagging. The tag integrated within steri-ID instrument baskets is not only compatible with various suppliers of RFID writing and reading technologies, but can be scanned through sterilisation packaging. The integrated tag is available in two standard sizes (DIN 1/1 and DIN 1/2).

The Steri-ID range is also resistant to all common washing and sterilisation conditions, supplied in 5 different colours to distinguish between distinct surgical sets and can be implemented for over 400 cycles.



CBM Baskets



The CBM basket ranges are manufactured with AISI304 grade stainless steel and can be introduced within sterilisation containers or with sterilisation barrier wrap. These baskets are also resistant to colour changes as a result of electro polishing and exhibit a surface treatment that increases their water drainage capacity.

Additional benefits include a tailored design service which enables baskets or silicone supports to be produced in line with user requirements. Please contact a member of our sales team for further size enquiries.



MINI AND MICRO BASKETS

CODE	DESCRIPTION	SIZE
CBM63J	Mini basket with lid (without handles)	200 x 80 x 35mm
CBM63M	Mini basket with lid (without handles)	220 x 50 x 35mm
CBM67F2	Micro basket with lid (without handles)	170 x 135 x 35mm
CBM70F2	Micro basket with lid and handles	210 x 145 x 50mm
CBM76D	Mini basket with lid and handles	267 x 125 x 49mm
CBM77F2	Micro basket with lid and handles	270 x 170 x 35mm



SCOPE BASKETS SYSTEM

CODE	DESCRIPTION	SIZE
CBM48/S	Scope basket for 1 scope (with lid)	480 x 150 x 80mm
CBM66/S	Scope basket for 1 scope (with lid)	450 x 70 x 60mm
CBM65/SA	Scope basket for 1 scope (with lid)	290 x 70 x 60mm
CBM64ECV	Empty scope basket (with lid and handles)	650 x 150 x 70mm
CBM66AV	Empty scope basket (with lid and handles)	580 x 150 x 80mm



STERILISATION BASKETS SYSTEM

CODE	DESCRIPTION	SIZE
CBM68/1	Stainless steel basket (without lid)	340 x 240 x 50mm
CBM68LID	Stainless steel basket (with lid)	340 x 240 x 50mm
CBM71/1	Stainless steel basket (without lid)	240 x 240 x 50mm
CBM71LID	Stainless steel basket (with lid)	240 x 240 x 50mm
CBM72/1	Stainless steel basket (without lid)	300 x 240 x 35mm
CBM72LID	Stainless steel basket (with lid)	300 x 240 x 35mm
CBM74/1	Stainless steel basket (without lid)	480 x 240 x 50mm
CBM74LID	Stainless steel basket (with lid)	480 x 240 x 50mm
CBM75/1	Stainless steel basket (without lid)	540 x 240 x 50mm
CBM75LID	Stainless steel basket (with lid)	540 x 240 x 50mm
CBM79/2	Stainless steel basket (without lid)	650 x 240 x 70mm
CBM79LID	Stainless steel basket (with lid)	650 x 240 x 70mm

Plastic scope trays



Plastic scope trays facilitate the sterilisation, organisation and storage of delicate instruments. Their material will not damage instruments and is compatible with all standard sterilisation methods.



CODE	PRODUCT	DESCRIPTION	SIZE
P261615	Scope tray	Base, lid & inserts	66 x 406 x 38mm
P651815	Scope tray	Base, lid & inserts	165 x 457 x 38mm
P651815M	Scope tray	Base, lid & mats	165 x 457 x 38mm
P651830P3	Scope tray	Base, lid & inserts	165 x 457 x 76mm
P651830/HED	Scope tray	Base, lid & mats	165 x 457 x 76mm
P6520A	Scope tray	Base, lid & inserts	76 x 559 x 38mm

Stainless steel scope baskets



A lockable stainless steel wire mesh basket for single endoscopes/arthroscopes with integral silicone instrument holders. This product range protects expensive scope instruments from the rigours of washing and sterilising.



FOR ONE ENDOSCOPE

CODE	SIZE
CIL2728	290 x 80 x 52mm
CIL2727	460 x 80 x 52mm
CIL2726	670 x 80 x 52mm
CIL2701	640 x 150 x 77mm (with small mesh basket)

FOR TWO ENDOSCOPES

CODE	SIZE
CIL2730	290 x 80 x 52mm
CIL2729	460 x 80 x 52mm

Silicone inserts



Silicone inserts allow custom configuration within the tray to accept various sizes of instruments.

CODE	SIZE
A387592	95 x 19mm
A387593	95 x 19mm
A387594	95 x 19mm
A387595	95 x 19mm
A387596	95 x 19mm
A237591K	57 x 13mm
A237592K	57 x 13mm
A237593K	57 x 13mm
A237593	57 x 13mm




CODE	SIZE
A237596	57 x 13mm
A3875912	95 x 32mm
A32759P1	95 x 32mm
A32759P2	95 x 32mm
A327593	95 x 32mm
A327593L	95 x 32mm
A2616901	53 x 25mm
A2616902	53 x 25mm

Delicate wash tray baskets



Delicate wash tray baskets are manufactured with stainless steel and silicone rests have been fitted into their base. Instruments sit in individual rests and are held in place with a firm clip protected with a silicone shield. Fine instruments are particularly protected during washing and sterilisation.



CODE	PRODUCT	DESCRIPTION	SIZE
RON005	Full tray	Can hold 20 instruments	285 x 185mm
RON006	Half tray	Can hold 9 instruments	142 x 185mm
RON007	Quarter tray	Can hold 5 instruments	71 x 185mm

Aluminium trays



A range of trays which is manufactured with 16swg aluminium sheets with folded edges for extra rigidity and boasts welded or open corners with side holes to ensure tray tags can be easily fixed. Lids are available to fit trays.






CODE	SIZE
TRY001	305 x 152 x 51mm
TRY002	305 x 305 x 51mm
TRY003	305 x 457 x 51mm
TRY004	305 x 610 x 51mm
TRY005	457 x 457 x 51mm
TRY201	356 x 495 x 51mm (open cornered)

Plastic sterilisation trays



A range of trays designed to cover most requirements for protecting instruments during sterilisation and subsequent transportation. The trays exhibit grid holes to provide drainage, fixing positions for various inserts, and buttons to raise silicone pin mats from the base. Transparent lids are strongly fitted to provide additional security and ensure all contents are highly visible. Steam, ethylene oxide or formaldehyde sterilisation processes are all suitable.



CODE	PRODUCT	DESCRIPTION	SIZE
P0278	Mini tray	Mini tray (base, lid & mat)	38 x 68 x 19mm
P0278125	Mini tray	Mini tray (base, lid & mat)	38 x 68 x 32mm
P0625	Microsurgical instrument tray	Micro tray (base, lid & mat)	64 x 125 x 19mm
P0625K1	Microsurgical instrument tray	Micro tray (base, lid & mat) – 2 silicone inserts with 1 instrument slot	64 x 152 x 19mm
P0625K3	Microsurgical instrument tray	Micro tray (base, lid & mat) – 2 silicone inserts with 3 instrument slots	64 x 152 x 19mm
P2610A	Microsurgical instrument tray	Micro tray (base, lid & mat)	64 x 125 x 32mm
P0725	Microsurgical instrument tray	Micro tray (base, lid & mat)	64 x 190 x 19mm
P6540	Microsurgical instrument tray	Micro tray (base, lid & mat)	102 x 162 x 19mm
P7540	Microsurgical instrument tray	Micro tray (base, lid & mat)	102 x 190 x 19mm
P100675	Microsurgical instrument tray	Micro tray (base, lid & mat)	152 x 254 x 19mm
P1006151	Microsurgical instrument tray	Micro tray (base, lid & mat)	152 x 254 x 38mm
P1006153	Microsurgical instrument tray	Micro tray (base, lid & mat)	152 x 254 x 19mm
P101510	Microsurgical instrument tray	Micro tray (base, lid & mat)	254 x 381 x 25mm
P1015251	Microsurgical instrument tray	Micro tray (base, lid & mat)	254 x 381 x 25mm

Healthmark trays



A selection of heat and chemical resistant polypropylene trays from the USA with perforations that enable silicone inserts to be effortlessly introduced. These trays are perfect for lower budgets and most ideal when all inner contents do not need to be viewed throughout the sterilisation process.



CODE	PRODUCT	SIZE
HM30732-H	Healthmark tan economy tray	184 x 83 x 29mm
HM30842-H	Healthmark tan economy tray	222 x 114 x 44mm
HM31042-H	Healthmark tan economy tray	260 x 114 x 44mm
HM31462-H	Healthmark tan economy tray	356 x 159 x 64mm
HM31772-H	Healthmark blue economy tray	445 x 197 x 60mm



Inspection

pages 17-22

Absorbent towels

Low-linting towels designed to remove excess moisture from trays/instruments before the packing procedure.

Veraclean low-linting absorbent towels (Previously known as Chux)



Low-linting super absorbent single use crepe towels which can be used in the wash room to remove moisture from manually cleaned instruments and medical devices. These towels can also be used within the inspection, assembly and packing room to remove excess moisture droplets from instruments following automated cleaning and disinfection.



CODE	PRODUCT	QUANTITY
CHUX-WIPE01	Chux absorbent wipes (34 x 30cm)	400
CHUX-DISP	Chux absorbent wipe dispenser	1

Identification

A range of innovative products to not only identify instruments that require additional actions (such as repairs or sharpening), but to highlight essential instructions in a CSSD unit.

Instrument locks



A simple range of special locks designed for marking instruments and trays. The locks are printed with 'good', 'repair' and 'sharpen'. Specifically for use with instruments, for circulating between the user, the CSSD and the repair shops.



CODE	TEXT ON TAG
LOCK994G	Good
LOCK994R	Repair
LOCK994S	Sharpen

Washer tags



A range of self-tie tags which can be fixed onto the handles of wash baskets to communicate information to the clean room. The material can be written on with suitable marker pen or supplied pre-printed (subject to special order).



CODE	DESCRIPTION	SIZE	QUANTITY
TAG-WASH	Washer tags (white)	50 x 228.6mm	500 per roll
TAG-WASH-ORA	Washer tags (orange)	50 x 228.6mm	500 per roll
TAG-WASH-YEL	Washer tags (yellow)	50 x 228.6mm	500 per roll
TAG-WASH-GRN	Washer tags (green)	50 x 228.6mm	500 per roll
TAG-WASH-RED	Washer tags (red)	50 x 228.6mm	500 per roll
TAG-WASH-BLU	Washer tags (blue)	50 x 228.6mm	500 per roll
TAG-WASH-PNK	Washer tags (pink)	50 x 228.6mm	500 per roll
TAG-CLEANUSE	Washer tag printed with the words 'Clean ready for use'	50 x 228.6mm	500 per roll

Tracking & tracing

A selection of tracking and tracing controls to facilitate the identification of specific containers, trays and instruments. All controls consist primarily in colour coding, attachments and adhesive labelling.

Tray tags



Clinipak tray tags not only simplify the identification of sets or instruments, but provide a useful means of labelling items for tracking and tracing. All tags are produced with a polypropylene polymer that incorporates glass nylon to withstand the rigours of the decontamination cycle. Available in 17 colours and also in stainless steel.



CODE	COLOUR	SIZE	QUANTITY
TAG011	White	90 x 38mm	100
TAG012	Tan	90 x 38mm	100
TAG013	Yellow	90 x 38mm	100
TAG014	Orange	90 x 38mm	100
TAG015	Red	90 x 38mm	100
TAG016	Blue	90 x 38mm	100
TAG017	Green	90 x 38mm	100
TAG018	Grey	90 x 38mm	100
TAG019	Purple	90 x 38mm	100
TAG020	Black	90 x 38mm	100
TAG021	Lime	90 x 38mm	100
TAG022	Pink	90 x 38mm	100
TAG023	Turquoise	90 x 38mm	100
TAG024	Mustard	90 x 38mm	100
TAG025	Terracotta	90 x 38mm	100
TAG026	Royal blue	90 x 38mm	100
TAG027	Sky blue	90 x 38mm	100
TAG55WH	White	90 x 55mm	100
TAGMET7530	Stainless steel	75 x 30mm	20

Security ties for tray tags



CODE	DESCRIPTION	SIZE	QUANTITY
TIE001	Blue tie	3.6 x 144mm	100
TIE002	Clear tie	2.5 x 150mm	100

Split rings



CODE	SIZE	QUANTITY
RIN001	24mm stainless steel	100

Marking tape



Clinipak marking tape ensures that specified instruments can be returned to the correct department or surgeon after sterilisation and eradicates the requirement for engraving. The colour strips are long-lasting and appropriate for steam sterilisation. Independent research confirms that no bacteria remain under the tapes after sterilisation.



RING BINDERS

CODE	DESCRIPTION
MKT001	Ring binder with index only
MKT002	Ring binder with a sheet of plain colours (11 colours)
MKT003	Ring binder with a sheet of zebra colours (10 colours)
MKT004	Ring binder with a one sheet of each plain and zebra colours (21 colours)

TIGER TAPES (colour with black stripes)

CODE	DESCRIPTION	QUANTITY
MKT030	Yellow	5 sheets per pack
MKT031	Orange	5 sheets per pack
MKT032	Red	5 sheets per pack
MKT033	Green	5 sheets per pack

PLAIN TAPES

CODE	DESCRIPTION	QUANTITY
MKT005	Purple	5 sheets per pack
MKT006	White	5 sheets per pack
MKT007	Tan	5 sheets per pack
MKT008	Yellow	5 sheets per pack
MKT009	Orange	5 sheets per pack
MKT010	Red	5 sheets per pack
MKT011	Blue	5 sheets per pack
MKT012	Green	5 sheets per pack
MKT013	Silver	5 sheets per pack
MKT014	Brown	5 sheets per pack
MKT015	Black	5 sheets per pack

ZEBRA TAPES (colour with white stripes)

CODE	DESCRIPTION	QUANTITY
MKT016	Beige/Tan	5 sheets per pack
MKT017	Yellow	5 sheets per pack
MKT018	Orange	5 sheets per pack
MKT019	Red	5 sheets per pack
MKT020	Blue	5 sheets per pack
MKT021	Green	5 sheets per pack
MKT022	Silver	5 sheets per pack
MKT023	Brown	5 sheets per pack
MKT024	Blue	5 sheets per pack
MKT025	Purple	5 sheets per pack

Identification tags

S HP W 2 LF

Clinipak identification tags comply with different tracking and tracing systems and sterilisation methods. All tags and labels are supplied in plain white as standard, but are also available in different colours or pre-printed at an additional cost. The following colours are available: blue, green, yellow, red and orange.

CODE	PRODUCT	DESCRIPTION	SIZE	QUANTITY
SELFTIE01	Self-tie	Plain tag with dumbbell self-tie. Laserplas 195 mic.	25 x 190.5mm	1000 per roll
TAG-CONTAIN1	Small container tag	Self-tie tag with pear-drop tie. HTP 175mic.	75 x 228.6mm	250 per roll
TAG-CONTAIN2	Large container tag	Self-tie tag with pear drop tie. HTP 175mic.	95 x 279.4mm	250 per roll
TAG-SECURE1	Security tag	Plain white container security tag. 175 mic.	101.6 x 152mm	500 per roll
TAG-TRAY1	Tray tag	Plain tag which can be attached to a basket or tray using a cable tie. Laserplas 127mic.	70 x 101.6mm	500 per roll
TAG-AUTO1	Autoclave tag	Plain white tag with a t-shaped tie.	50 x 228.6mm	500 per roll



TAG-CONTAIN1



TAG-CONTAIN2



TAG-SECURE1



TAG-TRAY1



TAG-AUTO1



SELFTIE01

Marker pens

LF

Clinipak permanent marker pens are effective for writing on the surface of Clinipak tray tags, pouches and reels as their ink is autoclavable, non-toxic and quick-drying.



CODE	DESCRIPTION
PEN001	Black marker pen
PEN002	Red marker pen
PEN003	Blue marker pen

Numbered instrument locks

S EO HP ~~2~~ LF

Clinipak numbered instrument locks display unique numbers to aid the tracking, tracing and record-keeping processes. These locks are available in green (with or without numbers), red and yellow.



CODE	DESCRIPTION	QUANTITY
SNAPLOCKG001	Green (with numbers)	100
SNAPLOCKG002	Green (without numbers)	100
SNAPLOCKR002	Red (with numbers)	100
SNAPLOCKY002	Yellow (with numbers)	100

Self-adhesive Alert Labels

S W ~~2~~

A growing range of self-adhesive labels which have been designed to highlight particular requirements and facilitate the identification of wrapped sets in a sterile services department. All Clinipak labels are autoclaveable, demonstrate strong adhesion, utilise bright colours where possible and are applied to various types of outer sterilisation wrap.



CODE	DESCRIPTION	SIZE	QUANTITY
SAL-FASTTRACK	Semi-gloss pre-printed and self-adhesive red label with white 'FAST TRACK' lettering to indicate urgency.	76 x 76mm	1000 labels per roll
SAL-HEAVYSET	Semi-gloss pre-printed and self-adhesive yellow label with black 'HEAVY SET' lettering to indicate the need for careful lifting.	76 x 76mm	1000 labels per roll
SAL-LOANSET	Semi-gloss pre-printed and self-adhesive green label with white 'LOAN SET' lettering to indicate that a particular set is on loan.	76 x 76mm	1000 labels per roll
SAL-FRAGILE	Semi-gloss pre-printed and self-adhesive white label with red 'FRAGILE INSTRUMENTS' lettering to urge careful handling.	76 x 76mm	1000 labels per roll
SAL-DONOTSTACK	Semi-gloss pre-printed and self-adhesive white label with red 'DO NOT STACK' lettering.	76 x 76mm	1000 labels per roll
SAL-DISINFECT	Semi-gloss pre-printed and self-adhesive white label with red 'DISINFECTED ONLY DO NOT AUTOCLAVE' lettering.	76 x 76mm	1000 labels per roll



Packing

pages 23-38

Our flexible production methods ensure the best sizes can be produced for our customers, even if these are not the industry standard. All sizes can be manufactured smaller without sacrificing high performance or protection levels. Clinipak actively engages with customers to identify when sizes can be reduced, which results in cost savings, lower disposal costs and diminished waste.

Some of our most popular size options are in the table below, but please contact a member of our sales team for further size enquiries.

50 x 50cm	90 x 90cm	100 x 120cm	120 x 150cm	150 x 150cm
60 x 60cm	100 x 75cm	100 x 140cm	130 x 150cm	170 x 150cm
75 x 75cm	100 x 100cm	120 x 120cm	140 x 140 cm	180 x 150cm

Barrier wraps

A wide range of sterilisation wraps which provide a bacterial barrier and protect the contents of a set or pack from contamination. If classed as alcohol repellent, our barrier wraps can also be used within the sterile field.

All Clinipak barrier wraps conform to both ISO11607-1 and EN868-2.

Crepe and plain paper

A complete range of cellulose sterilisation wraps.

S EO G ~~2~~ CE LF Custom sizes available



- 60gsm
- ISO11607-1 and EN868-2 compliant
- Available in green or white

Clinipak 340 wrap

Clinipak 340 wrap is a high performance sterilisation material that is suitable for inner and outer wrapping of medium to large trays and packs. This material exhibits a synthetic reinforcement within the cellulose fibre structure, which improves strength and increases its softness.

S EO ~~2~~ CE LF Custom sizes available



- 60gsm
- ISO11607-1 and EN868-2 compliant
- Alcohol repellent
- Available in green and blue
- Also available for Gamma sterilisation

Our wrap range is also designed to overcome many of the common challenges encountered when sterilising, storing and transporting wrapped packs.

Choice wrap

S EO HP 2 CE LF Custom sizes available

A fully synthetic SMS product which has been designed to offer total flexibility and enable users to pack conventionally or with 'two layers' at once. Clinipak choice wrap is not just exceptionally strong and very drapeable, but is available with alcohol repellence in both green and blue for use within the sterile field. The blue wrap can be supplied without alcohol repellency if required.

A 70gsm choice blue ultra edition of this product is also available to be used as an outer sterilisation wrap, given that its high strength provides effective protection against cuts and tears, particularly against heavy orthopaedic sets.

All choice wrap editions comply with the essential requirements stated within the Medical Devices Directive EC 93/42 - 2007/47/EC - Class 1 and are CE-marked on their product labels.



- 53 – 55gsm
- ISO11607-1 and EN868-2 compliant
- Available in green or blue
- Green wrap is alcohol repellent
- Blue wrap is available with OR without alcohol repellency
- Also suitable for Hydrogen Peroxide/ Plasma sterilisation

Choice wrap dry

S EO 2 CE LF Custom sizes available

A patented tray wrap solution to eliminate or heavily reduce the occurrence of wet packs. Clinipak choice wrap dry differs from conventional wraps, given that the inner surface exhibits a high degree of absorbency and the outer surface repels both water and alcohol. The absorbent surface disperses the condensate within a typical steam sterilisation process and assists with evaporation during the vacuum cycle.

This wrap is manufactured with a blend of materials that are compliant with the recognised standards for toxicology, bacterial resistance and drape.



- 79gsm
- ISO11607-1 and EN868-2 compliant
- Alcohol repellent
- Available in green or blue

Wrap bags

Both barrier and transportation wraps are available in a 'pillowcase' format to facilitate the wrapping process. Simply fold and seal the bags with autoclave tape for a time-saving alternative. Custom sizes are available to suit distinct requirements.

Choice bags

An innovative solution that enables users to introduce a strong, autoclavable outer bag as a transportation bag or a second layer of sterilisation wrap. Manufactured from the same strong material as choice wrap, Clinipak choice bags provide an outer layer that is incredibly user-friendly. Wrapped trays can be easily inserted into the bag, closed and sealed with autoclave tape.

S **EO** **HP** ~~2~~ **CE** **LF** Custom sizes available

- 55gsm
- ISO11607-1 and EN868-2 compliant
- Available in green or blue
- Green wrap is alcohol repellent
- Blue wrap is available with OR without alcohol repellency
- Also suitable for Hydrogen Peroxide/Plasma sterilisation

Choice dry bags

Manufactured from the same material as choice wrap dry, choice dry bags assist the sterilisation process in the same way to enable the removal of condensate and encourage dry packs. Wrapped trays can be easily inserted into the bag, closed and sealed with autoclave tape.

S **EO** ~~2~~ **CE** **LF** Custom sizes available



CODE	COLOUR	SIZE	QUANTITY
CDBAG001B	Blue	400 x 700 x 90mm (Lip)	200
CDBAG001G	Green	400 x 700 x 90mm (Lip)	200
CDBAG002B	Blue	500 x 700 x 90mm (Lip)	200
CDBAG002G	Green	500 x 700 x 90mm (Lip)	200
CDBAG003B	Blue	600 x 700 x 90mm (Lip)	100
CDBAG003G	Green	600 x 700 x 90mm (Lip)	100

Transportation cover bags

All transportation wraps are also available in a bag format in which wrapped trays can be effortlessly inserted. Clinipak transportation cover bags are closed using the provided lip and fastened with autoclave tape.

S **EO** ~~2~~ **CE** **LF** Custom sizes available



CODE	COLOUR	SIZE	QUANTITY
TCOVER001	Purple	400 x 700 x 90mm (Lip)	200
TCOVER002	Purple	500 x 700 x 90mm (Lip)	200
TCOVER003	Purple	600 x 700 x 90mm (Lip)	100
TCOVERDRY001	Purple Dry	400 x 700 x 90mm (Lip)	200
TCOVERDRY002	Purple Dry	500 x 700 x 90mm (Lip)	200
TCOVERDRY003	Purple Dry	600 x 700 x 90mm (Lip)	100

Transport wraps

Transportation wraps have been made to provide a strong dust cover and resist tearing. All Clinipak transportation wraps must be used in conjunction with one compliant barrier wrap to completely prevent the ingress of micro-organisms.

Use of a single barrier layer with transportation wrap can be a cost-effective solution which reduces damage in handling and transport. If adopted, changes to both wrapping and unwrapping procedures may be required to ensure aseptic opening.

Transportation cover wrap

S EO HP ~~Z~~ CE LF Custom sizes available

The ultimate handling material used to combat damaged packs due to its strength and tear resistance. The range is not only convenient when handling heavy orthopaedic sets, but uses distinctive colours to easily distinguish them from barrier layers. All Clinipak transportation wraps can be supplied on their own, interleaved with any compliant barrier layer or introduced as a triple layer to guarantee a strong, high-quality and cost-effective wrap combination.



- 60gsm
- Must be used in combination with a compliant barrier layer
- The range of colours can be used to differentiate packs e.g. loan sets, offsite decontamination etc.
- Available in purple, orange, yellow and red
- Also suitable for Hydrogen Peroxide/Plasma sterilisation

Dry transportation cover wraps

S EO ~~Z~~ CE LF Custom sizes available

Dry transportation cover wrap consists of purple transportation wrap with an additional, patented absorbent layer. This wrap provides optimum protection from dust, moisture and handling in addition to an absorbent inner layer to facilitate the dispersion of condensate during the autoclaving process.



- 75 and 83gsm
- Must be used in combination with the compliant barrier layer
- Can be used to identify specific packs
- Available in purple

Autoclave tapes

A wide range of unique, manufactured tapes which demonstrate excellent binding characteristics and comply with various types of sterilisation wraps.

Standard autoclave tapes

CODE	PRODUCT	DESCRIPTION	QUANTITY	SIZE	PROPERTIES
ST100.25NUFF	Steam autoclave tape with indicator	This has an effective indicator turning from white to dark brown/black when subject to the specified conditions.	36/ctn	25mm x 50m	S 2
ST130PL25	Plain steam autoclave tape	The same quality as the autoclave tape, but without an indicator. This tape is for closing purposes only.	36/ctn	25mm x 50m	S 2
ST101SMS19	SMS tape	A blue, steam autoclave tape which exhibits excellent adhesive properties for synthetic wraps etc.	48/ctn	19mm x 50m	S 2
ST.110ETO	Gas ethylene oxide tape	A tape that changes colour when exposed to ethylene oxide gas.	48/ctn	19mm x 50m	EO 2

Latex-free tapes

CODE	PRODUCT	DESCRIPTION	QUANTITY	SIZE	PROPERTIES
LF19TAPE3311	Regular latex free tape 19mm	This has the same properties as our steam autoclave tape, but latex free.	48/ctn	19mm x 50m	S 2 LF
LF25TAPE3311	Regular latex free tape 25mm	This has the same properties as our steam autoclave tape, but latex free.	36/ctn	25mm x 50m	S 2 LF
LFWPTAPE19	Heavy duty latex free tape 19mm	This tape has the same properties as our steam autoclave tape, but benefits from stronger adhesion. This latex-free tape is often used with synthetic-type wraps.	48/ctn	19mm x 50m	S 2 LF
LFWPTAPE25	Heavy duty latex free tape 25mm	This tape has the same properties as our steam autoclave tape, but benefits from stronger adhesion. This latex-free tape is often used with synthetic-type wraps.	36/ctn	25mm x 50m	S 2 LF
LF25TAPE3390PLAS	Plasma tape	This tape is suitable to be used in plasma sterilisation processes.	48/ctn	25mm x 25m	HP 2 LF

Tape dispensers

A heavy duty tape dispenser for use on workbenches.

CODE	PRODUCT
DISP004	3 reel heavy duty dispenser
DISP005	2 reel heavy duty dispenser




Protective packaging

All Clinipak packing products are designed to protect high-cost equipment from damage during the autoclave process and in subsequent transport and handling. These products are suitable for steam sterilisation, form a strong layer for additional safety, and can be effortlessly introduced at the time of packing.

Clinipak Corners

A simple and effective solution to reduce damage from tray corners and feet.

CODE	PRODUCT	DESCRIPTION	QUANTITY	PROPERTIES
CC001	Card corners	Clinipak corners rigid card	200/ctn	S EO 2 CE LF
CC002	Crepe corners	Clinipak corners absorbent crepe	100/ctn	S EO 2 CE LF
STOP-7100	Silicone corners	Green silicone corner and edge protector	40/ctn	S EO CE LF
STOP-7100P500	Silicone corners	Green silicone corner and edge protector	500/ctn	S EO CE LF



Tray anchors

Clinipak tray anchors provide optimum protection as instruments are tightly locked into place within the basket or tray after washing. This ensures that instruments do not move during subsequent handling and the potential for damage is significantly reduced. The material can be cut to length as required.



CODE	QUANTITY	PROPERTIES
ANC001	16	S 2 CE LF
ANC002	64	S 2 CE LF

Point protectors



Clinipak point protectors protect sharp instruments and needles during sterilisation, transport and storage. Effective, economical and safe to discard. 'Scissor' type instruments may be inserted, opened and locked open by folding the top cut upwards to ensure effective steam circulation. Three varieties are available in two different pack quantities: needle, 15mm & 30mm width slits.



POINT PROTECTORS

CODE	SIZE	QUANTITY
PPT001	15mm wide	2000
PPT002	15mm wide	1000
PPT003	30mm wide	2000
PPT004	30mm wide	1000

NEEDLE PROTECTORS

CODE	QUANTITY
PPT006	1000

Transparent point protectors



Transparent point protectors protect the tips and points of specific instruments and facilitate their visual identification. Their material is fully autoclavable and suitable for steam, plasma and ethylene oxide sterilisation. The unique design minimises the movement of instruments and holds them open to strengthen their subsequent sterilisation.



CODE	SIZE	QUANTITY
PPT007C	15mm wide	200
PPT008C	30mm wide	200

Chessboard instrument protector



An autoclave foam sheet with chessboard design cuts. All pieces can be removed to create individual nests for each instrument. The chessboard instrument protector is supplied with a top cover sheet.



CODE	SIZE
FOB003	305 x 457 x 55mm

Foam instrument pouches



Clinipak instrument pouches are produced with a high-quality autoclavable soft foam to ensure all instruments are delicately protected and grouped together. Each pouch has 12 sections.



CODE	QUANTITY
FOP001	1 single pouch
FOP002	50 x FOP001

Silicone pin matting



A silicone pin mat with drainage holes and buttons on the underside to allow for easy drainage and steam penetration. The material can withstand all known washing methods and can be cut to size as required.



CODE	SIZE
DINPIN001	500 x 300mm

Paper multipouches



Paper multipouches are designed to provide protection for individual instruments during sterilisation and transport. This product is manufactured from paper to ensure steam penetration when used in containers.



CODE	SIZE	QUANTITY
WMPP40X130	40 x 130mm	1000
WMPP77X130	77 x 130mm	1000

Absorbent liners

An extensive range of liners designed to absorb and spread condensate during autoclaving to encourage drier packs. As a result of our flexible production methods, bespoke sizes can be manufactured on request.

Millennium liners (40gsm)

S EO 2 CE LF Custom sizes available

Clinipak millennium blue liners are lightweight, user-friendly and produced with a matt blue finish to minimise glare. The liners are a 100% viscose-based dry laid product with a similar absorbency per unit area as our tray liners (80gsm).

Clinipak can manufacture almost any size up to a maximum width of 500mm x any length, but a selection of our most popular sizes can be seen in the table below.



CODE	SIZE	QUANTITY
MIL002	310 x 410mm	1000
MIL004	240 x 250mm	2400
MIL005	250 x 500mm	1200
MIL007	320 x 250mm	2000
MIL009	350 x 500mm	1200
MIL011	320 x 460mm	1000
MIL016	450 x 610mm	1000
MIL019	305 x 305mm	2000
MIL020	500 x 500mm	1000
MIL022	305 x 610mm	1000

Tray liners (80gsm)

S EO 2 CE LF Custom sizes available

Clinipak tray liners (80gsm) liners are produced with a finely creped medical quality paper. The high absorbency and rapid dispersion ensures condensate is removed during autoclaving, which encourages dry packs and improves the effectiveness of the sterilisation process. Tray liners can be used above or beneath all trays to protect against sharp basket corners.

Clinipak can manufacture almost any size up to a maximum width of 457mm x any length, but a selection of our most popular sizes can be seen in the table below.



CODE	SIZE	QUANTITY
TLC005	250 x 250mm	2400
TLC006	250 x 500mm	1200
TLC008	457 x 550mm	800
TLC009	305 x 305mm	2000
TLC010	305 x 610mm	1000
TLC011	305 x 457mm	1000
TLC013	457 x 610mm	600
TLC055	250 x 550mm	1200
TLC085	250 x 480mm	1000
TLC110	150 x 250mm	2000

Cut corner liners

S EO 2 CE LF Custom sizes available

80gsm crepe tray liners with pre-cut corners to improve the process of fitting the liner to the tray. These tray liners cover all surfaces of the tray including the sides.

Soaker sheets (140gsm)

S EO 2 CE LF Custom sizes available

Clinipak soaker sheets are used under trays to not only protect inner wraps from the sharp corners or feet on trays and baskets, but to absorb condensate during sterilisation. These sheets are manufactured from a heavier medical paper than our tray liners and can be used as an shelf liner for added protection from rough handling and wet packs.

Clinipak can manufacture almost any size up to a maximum width of 750mm x any length, but a selection of our most popular sizes can be seen in the table below.



CODE	SIZE	QUANTITY
SOK001SS	750 x 500mm	50
SOK002SS	500 x 375mm	100
SOK003SS	305 x 610mm	100
SOK004SS	250 x 500mm	100
SOK005SS	500 x 500mm	100
SOK006SS	250 x 250mm	100
SOK007SS	280 x 500mm	100
SOK008SS	560 x 1219mm	50
SOK009SS	560 x 1524mm	50
SOK011SS	584 x 914mm	50

Absorbent and protective mats

S EO 2 CE LF Custom sizes available

Clinipak absorbent mats consist of a strong polypropylene sheet that can be placed under trays to protect against sharp edges. The 'dry' version, with an added absorbent layer, is also available to provide rapid moisture dispersion during steam sterilisation. The purple mats are suitable for hydrogen peroxide sterilisation processes and the green mats exhibit higher levels of absorbency.



PURPLE	PURPLE DRY	GREEN DRY	SIZE	QUANTITY
MAT001	DRYMAT001	DRYMATGRN001	750 x 500mm	50
MAT002	DRYMAT002	DRYMATGRN002	500 x 375mm	100
MAT003	DRYMAT003	DRYMATGRN003	305 x 610mm	100
MAT004	DRYMAT004	On request	250 x 500mm	100
MAT005	DRYMAT005	On request	250 x 250mm	100
MAT006	DRYMAT006	On request	280 x 500mm	100
MAT007	DRYMAT007	On request	500 x 500mm	100
MAT008	DRYMAT008	On request	600 x 600mm	100
MAT009	DRYMAT009	On request	750 x 750mm	100
MAT010	DRYMAT010	On request	900 x 900mm	100

Bags, pouches & reels

A diverse range of sterilisation bags, pouches and reels in a wide variety of sizes to simplify the packing of surgical instruments. Clinipak dispensers ensure all products are organised and easy to locate.

Clinipak standard pouches & reels



Clinipak pouches and reels are produced with a 68gsm high burst resistance medical grade paper, which is sealed to tear resistant film to provide excellent seal strength and good peel characteristics. All standard items are suitable for steam and ETO sterilisation, but a further range is available for hydrogen peroxide sterilisation procedures. All print and process indicators are located outside of the inner pouch or reel area to avoid contact with contents.

Clinipak can manufacture pouches and reels with various materials for all common sterilisation methods and in extra sizes up to a maximum of 900mm x 3000mm.

FLAT POUCHES

CODE	PRODUCT	QUANTITY	SIZE
CLFP075X270	Flat pouch	1000/ctn	75 x 270mm
CLFP100X150	Flat pouch	1000/ctn	100 x 150mm
CLFP100X250	Flat pouch	1000/ctn	100 x 250mm
CLFP100X270	Flat pouch	1000/ctn	100 x 270mm
CLFP100X300	Flat pouch	1000/ctn	100 x 300mm
CLFP100X350	Flat pouch	1000/ctn	100 x 350mm
CLFP100X400	Flat pouch	1000/ctn	100 x 400mm
CLFP150X300	Flat pouch	1000/ctn	150 x 300mm
CLFP150X350	Flat pouch	1000/ctn	150 x 350mm
CLFP150X380	Flat pouch	1000/ctn	150 x 380mm
CLFP200X400	Flat pouch	1000/ctn	200 x 400mm
CLFP225X420	Flat pouch	500/ctn	225 x 420mm
CLFP250X350	Flat pouch	1000/ctn	250 x 350mm
CLFP250X380	Flat pouch	1000/ctn	250 x 380mm
CLFP250X485	Flat pouch	500/ctn	250 x 485mm
CLFP300X450	Flat pouch	500/ctn	300 x 450mm
CLFP300X500	Flat pouch	500/ctn	300 x 500mm
CLFP400X500	Flat pouch	500/ctn	400 x 500mm



FLAT REELS

CODE	PRODUCT	QUANTITY	SIZE
CLFR050	Flat reel	5/ctn	50mm x 200m
CLFR075	Flat reel	4/ctn	75mm x 200m
CLFR100	Flat reel	3/ctn	100mm x 200m
CLFR150	Flat reel	2/ctn	150mm x 200m
CLFR200	Flat reel	1/ctn	200mm x 200m
CLFR250	Flat reel	1/ctn	250mm x 200m
CLFR300	Flat reel	1/ctn	300mm x 200m
CLFR400	Flat reel	1/ctn	400mm x 200m



GUSSETED POUCHES

CODE	PRODUCT	QUANTITY	SIZE
CLGP100X300	Gusseted pouch	500/ctn	100mm x 50mm x 300mm
CLGP100X360	Gusseted pouch	500/ctn	100mm x 50mm x 360mm
CLGP100X500	Gusseted pouch	500/ctn	100mm x 50mm x 500mm
CLGP150X400	Gusseted pouch	500/ctn	150mm x 50mm x 400mm
CLGP150X460	Gusseted pouch	500/ctn	150mm x 50mm x 460mm
CLGP150X600	Gusseted pouch	500/ctn	150mm x 50mm x 600mm
CLGP200X400	Gusseted pouch	500/ctn	200mm x 50mm x 400mm
CLGP250X480	Gusseted pouch	250/ctn	250mm x 65mm x 480mm

GUSSETED REELS

CODE	PRODUCT	QUANTITY	SIZE
CLGR075	Gusseted reel	4/ctn	75mm x 38mm x 100m
CLGR100	Gusseted reel	3/ctn	100mm x 50mm x 100m
CLGR150	Gusseted reel	2/ctn	150mm x 50mm x 100m
CLGR200	Gusseted reel	1/ctn	200mm x 50mm x 100m
CLGR250	Gusseted reel	1/ctn	250mm x 65mm x 100m
CLGR300	Gusseted reel	1/ctn	300mm x 65mm x 100m

SELF-SEAL POUCHES

CODE	PRODUCT	QUANTITY	SIZE
CLSS090X200	Self-seal pouch	2000/ctn	90 x 200mm
CLSS090X230	Self-seal pouch	2000/ctn	90 x 230mm
CLSS130X385	Self-seal pouch	1000/ctn	130 x 385mm
CLSS135X255	Self-seal pouch	2000/ctn	135 x 255mm
CLSS190X330	Self-seal pouch	1000/ctn	190 x 330mm
CLSS200X385	Self-seal pouch	1000/ctn	200 x 385mm
CLSS250X440	Self-seal pouch	1000/ctn	250 x 440mm
CLSS300X440	Self-seal pouch	1000/ctn	300 x 440mm
CLSS305X381	Self-seal pouch	1000/ctn	305 x 381mm

Tyvek® pouches and reels

A range of Tyvek® pouches and reels are suitable for both ETO and Hydrogen Peroxide sterilisation methods.

DuPont™ and Tyvek® are DuPont trademarks or registered trademarks.

Steriking premium pouches & reels



Steriking pouches and reels possess a unique nine-layer film with seven polypropylene (PP) layers which guarantee high elasticity, strength and durability. All Steriking medical pouches also exhibit strong, factory-made seals which can withstand vigorous processing and handling. These pouches and reels are produced in line with clean-room condition requirements, use various medical grade films or breathable webs, and suit all common sterilisation methods.

The premium Steriking range supplied through Clinipak includes: flat pouches, flat reels, gusseted pouches, gusseted reels, and self-seal pouches. The Tyvek® variant of these pouches and reels is also available for plasma and gamma sterilisation. **Please contact a member of our sales team for further size enquiries.**

STERIKING FLAT POUCHES

CODE	PRODUCT	NHSSC CODE	SIZE	QUANTITY
WI-S7	Flat pouch	ETE641	205 x 270mm	600/ctn
WI-S8	Flat pouch	ETE637	100 x 300mm	1200/ctn
WI-S9	Flat pouch	ETE639	100 x 400mm	1200/ctn
WI-S10	Flat pouch	ETE642	205 x 400mm	600/ctn
WI-S13	Flat pouch	ETE615	150 x 400mm	1200/ctn
WI-S15	Flat pouch	ETE640	150 x 200mm	1200/ctn
WI-S16	Flat pouch	ETE643	250 x 380mm	1200/ctn
WI-S18	Flat pouch	ETE636	75 x 300mm	4200/ctn
WI-S19	Flat pouch	ETE638	100 x 350mm	1200/ctn
WI-S20	Flat pouch	ETE616	150 x 300mm	1800/ctn
WI-S21	Flat pouch	ETE617	300 x 500mm	600/ctn
WI-S23	Flat pouch	ETE635	75 x 230mm	1800/ctn



STERIKING FLAT REELS

CODE	PRODUCT	NHSSC CODE	SIZE	QUANTITY
WI-R39	Flat reel	ETF632	50mm x 200m	5/ctn
WI-R40	Flat reel	ETF633	75mm x 200m	4/ctn
WI-R41	Flat reel	ETF634	100mm x 200m	3/ctn
WI-R42	Flat reel	ETF635	150mm x 200m	2/ctn
WI-R43	Flat reel	ETF644	200mm x 200m	1/ctn
WI-R44	Flat reel	ETF636	250mm x 200m	1/ctn
WI-R45	Flat reel	ETF637	300mm x 200m	1/ctn
WI-R46	Flat reel	ETF638	350mm x 200m	1/ctn
WI-R47	Flat reel	ETF645	400mm x 200m	1/ctn



STERIKING GUSSETED POUCHES

CODE	PRODUCT	NHSSC CODE	SIZE	QUANTITY
WI-B31	Gusseted pouch	ETE618	150mm x 50mm x 400mm	500/ctn
WI-B34	Gusseted pouch	ETE619	300mm x 80mm x 550mm	250/ctn
WI-B35	Gusseted pouch	ETE620	100mm x 50mm x 300mm	500/ctn
WI-B36	Gusseted pouch	ETE621	200mm x 55mm x 400mm	250/ctn
WI-B37	Gusseted pouch	ETE622	250mm x 65mm x 480mm	250/ctn

STERIKING GUSSETED REELS

CODE	PRODUCT	NHSSC CODE	SIZE	QUANTITY
WI-RB50	Gusseted reel	ETF639	75mm x 20mm x 100m	4/ctn
WI-RB51	Gusseted reel	ETF640	100mm x 50mm x 100m	3/ctn
WI-RB52	Gusseted reel	ETF646	150mm x 50mm x 100m	2/ctn
WI-RB53	Gusseted reel	ETF630	200mm x 55mm x 100m	1/ctn
WI-RB54	Gusseted reel	ETF641	250mm x 65mm x 100m	1/ctn
WI-RB55	Gusseted reel	ETF642	300mm x 80mm x 100m	1/ctn
WI-RB56	Gusseted reel	ETF643	350mm x 80mm x 100m	1/ctn
WI-RB57	Gusseted reel	ETF631	400mm x 80mm x 100m	1/ctn

STERIKING SELF-SEAL POUCHES

CODE	PRODUCT	NHSSC CODE	SIZE	QUANTITY
WI-SS1	Self-seal pouch	ETE623	90 x 200mm	1000/ctn
WI-SS2	Self-seal pouch	ETE624	90 x 250mm	1000/ctn
WI-SS4	Self-seal pouch	ETE625	130 x 270mm	1000/ctn
WI-SS4A	Self-seal pouch	ETE626	130 x 350mm	1000/ctn
WI-SS5	Self-seal pouch	ETE645	200 x 350mm	1200/ctn
WI-SS5A	Self-seal pouch	ETE644	190 x 330mm	1000/ctn
WI-SS6	Self-seal pouch	ETE627	250 x 400mm	600/ctn
WI-SS7	Self-seal pouch	ETE646	300 x 450mm	600/ctn

Dispensers & cutters

LF

A range of dispensers designed to organise access to bags, pouches or reels on a working surface.

CODE	PRODUCT
DISP001	Reel dispenser (supplied with cutter, pen holder and PEN001/PEN002/PEN003).
DISP002	Pouch and- bag rack (supplied with three large trays, one small tray, cutter, pen holder and PEN001/PEN002/PEN003).
CUT001	Hooked cutter to effectively cut flat and gusseted reels and open delivery packages.



* Reels and pouches not included

Sterilisation paper bags



A selection of plain and heat sealable paper sterilisation bags which has been produced with medical grade paper, adhesive and inks. This particular range includes a stepped top for easy filling, is available with or without heat seals, and has been effectively serrated to assist with health and safety. The sterilisation paper bags are compatible with all popular sterilisation methods and a plain top version is also available without indicators.

This product range is compliant with EN868 part 4.



CODE (Heat seal)	CODE (Plain top)	SIZE	QUANTITY
AHS	APC	90 x 50 x 150mm	1000
BHS	BPC	90 x 50 x 250mm	1000
CHS	CPC	110 x 30 x 190mm	1000
DHS	DPC	140 x 75 x 250mm	1000
EHS	EPC	190 x 65 x 250mm	1000
FHS	FPC	140 x 50 x 330mm	500
HHS	HPC	180 x 95 x 380mm	500
LHS	LPC	250 x 100 x 380mm	250
MHS	MPC	300 x 75 x 530mm	250
NHS	NPC	300 x 125 x 591mm	250
OHS	OPC	180 x 75 x 530mm	250
PHS	PPC	190 x 65 x 330mm	500
QHS	QPC	110 x 45 x 190mm	1000
RHS	RPC	125 x 50 x 250mm	250
T35HS	T35PC	380 x 125 x 610mm	250
THS	TPC	380 x 125 x 510mm	250
WHS	WPC	125 x 50 x 560mm	500



Sterilisation

pages 39-42

Cleaning and validation

Clinipak products for monitoring the performance of steam sterilisers and providing reliable means of validation.

Bowie & Dick packs



Bowie & Dick packs are single-use autoclave test packs to monitor the performance of steam sterilisers and provide a more reliable means of validation. This form of validation is fully compliant with European standard ISO 11140-4:2007. Easy to open and simple to implement. The test sheet **will not change** and consequently provides a safe and permanent record.



CODE	NHSSC CODE	QUANTITY
STE001	ETZ013	20

Container padlocks



Clinipak container padlocks are designed with a colour indicator that changes from blue to brown to indicate that the specified container has been processed by steam sterilisation. All padlocks can be easily fixed and have built-in break points to encourage easy removal. This product is compliant with ISO 11140-1.



CODE	DESCRIPTION	QUANTITY
LOCK3030BLUE	A small pack of container padlocks with a colour indicator.	100
LOCK3030BLUE1000	A large pack of container padlocks with a colour indicator.	1000

WIPAK multi-seal test kit



The WIPAK multi-seal test kit includes weekly and daily testing measures to monitor sealing processes. The testing consists of seal seam integrity testing, ISO 11607-1:2006, Striking® seal control sheets and seal leak detection with dye penetration ASTM F1929-12. This product enables users to test one sealing device for a one year period.



- 300 seal control test sheets
- 1 stopwatch
- 80ml dye solution (drop flacon and pipette)
- Instructions for use
- Weekly test reports
- Daily reports with sheets

WI SCAN

Coming soon!

WI SCAN is a patented, electronic Bowie & Dick test. This solution comprises of new air and NCG detection to analyse autoclave steam, which guarantees reliability in Bowie & Dick tests for all pre-vacuum autoclaves. WI SCAN simplifies controls, reduces the use of paper and optimises costs. **Ask a member of our sales team for more details.**

Sterilisation containers

A range of containers which protect and maintain the sterility of expensive instruments. All sterilisation containers ensure a clear return on investment and are a 'greener' alternative to sterilisation wraps.

CBM Containers



CBM Containers are installed with permaside aluminium filters, used in steam sterilizers conforming to EN285, and fully compliant with EN969-8 and ISO11607-1. The body is produced with anodized aluminium for increased corrosion resistance, particularly against modern alkaline washing cycles. The container lid is manufactured from stainless steel for improved scratch resistance. CBM containers also employ permanent PTFE filter cartridges which are certified for more than 2000 cycles, giving them a general shelf life of 5 years.




CODE	DESCRIPTION
CBM2070PS	600 x 300 x 260mm
CBM2071PS	600 x 300 x 210mm
CBM2072PS	600 x 300 x 160mm
CBM2073PS	600 x 300 x 110mm
CBM2079PS	600 x 300 x 135mm
CBM4072PS	460 x 300 x 160mm

CBM Container Filters



Clinipak distributes a specialised range of CBM container filters for steam sterilisation. In addition to single-use paper filters, CBM cotton filters and synthetic filters are available. The cotton filters can be implemented for over 70 cycles while the synthetic filters can be introduced for 200 cycles. **Please contact a member of our sales team for all size enquiries.**



polySteribox®



polySteribox® is a range of all-in-one, re-usable and hermetically-sealable sterilisation containers from Ritter Medical. The base demonstrates effective bacteria filters and the patented interlocking system prevents accidental opening. The containers effectively save space, can be easily stacked and are well-suited to both the storage and transportation of sterile goods. There are 2 forms of PTFE filter sets available: permanent/lifetime and annual.

Silicone mats and safety seals are also available exclusively for polySteribox®



CODE	DESCRIPTION
STERIBOX-SH	208 x 111 x 55mm
STERIBOX-M	185 x 145 x 59mm
STERIBOX-L	290 x 184 x 74mm
STERIBOX-XL	497 x 122 x 65mm

SUPPLEMENTARY PRODUCTS

CODE	PRODUCT
STERIBOX-MATSH	Silicone mat for polySteribox® SH
STERIBOX-SEAL	Safety seals for polySteribox® containers



Transport

pages 43-48

Transportation boxes and crates

A selection of boxes and crates which ensure that cleaned or contaminated surgical instruments can be transported effectively. These solutions reduce damage to wrapped sets in handling and transport.

Transportation boxes (UN3291)

Clinipak transportation boxes are suitable for the transportation of Medical Waste N.O.S (Class 6.2, Classification code I3, UN3291) in accordance with Packing Instruction P621 defined in ANNEX A of the European Agreement concerning the International Carriage of Dangerous Goods by Road. These transportation boxes are also compliant with section 2.2.62.1.5.7 of the agreement and facilitate the transportation of all used medical instruments or equipment for the purposes of decontamination, cleaning, sterilisation or repair.

These transportation boxes are stackable, can be custom printed with text on the box and lid, and withstand washer disinfectors for decontamination purposes. Blue and red variants of these boxes are available to distinguish between 'clean' and 'dirty' contents. Optional metal inner baskets are available for box sizes 002, 003 and 004.



BLUE TRANSPORTATION BOXES

W LF

CODE	DESCRIPTION	SIZE
BTB001	Weight: 0.6kg Capacity: 4 litres	External Size: 300 x 200 x 120mm Internal Size: 285 x 162 x 92mm
BTB002	Weight: 0.9kg Capacity: 7 litres	External Size: 300 x 200 x 170mm Internal Size: 287 x 162 x 153mm
BTB003	Weight: 1.4kg Capacity: 16 litres	External Size: 400 x 300 x 220mm Internal Size: 287 x 232 x 200mm
BTB004A	Weight: 3.3kg Capacity: 49 litres	External Size: 600 x 400 x 330mm Internal Size: 458 x 318 x 311mm
BTB005	Weight: 10.2kg Capacity: 170 litres	External Size: 800 x 600 x 600mm Internal Size: 568 x 489 x 570mm



RED TRANSPORTATION BOXES

W LF

CODE	DESCRIPTION	SIZE
RTB002	Weight: 0.9kg Capacity: 7 litres	External Size: 300 x 200 x 170mm Internal Size: 287 x 162 x 153mm
RTB003	Weight: 1.4kg Capacity: 16 litres	External Size: 400 x 300 x 220mm Internal Size: 287 x 232 x 200mm
RTB004	Weight: 3.3kg Capacity: 49 litres	External Size: 600 x 400 x 329mm Internal Size: 457 x 318 x 293mm
RTB006	Weight: 4.2 kg Capacity: 60 litres	External Size: 790 x 400 x 320mm Internal Size: 645 x 318 x 284mm



INNER BASKETS

S W LF

CODE	DESCRIPTION	SIZE
BTB002BASK	Inner basket to fit: BTB002/RTB002	155 x 275 x 145mm
BTB003BASK	Inner basket to fit: BTB003/RTB003	225 x 290 x 190mm
BTB004BASK	Inner basket to fit: BTB004/RTB004	440 x 310 x 285mm

Lidded crates (UN3291)

W LF

Clinipak lidded crates are a new addition to our range of boxes for the transportation of medical waste N.O.S. These crates include an integrated lid which can be sealed using either a pin or sealing strip to comply with transportation regulations. The range comprises of 4 sizes and boasts ergonomic handles for simple carrying. The two larger sizes (UNLC45 and UNLC60) can accommodate both wrapped and unwrapped 'DIN' sized trays



CODE	DESCRIPTION	SIZE
UNLC20	Weight: 1.68 kg Capacity: 20 litres	External Size: 400 x 300 x 275mm Internal Size: 325 x 227 x 238mm
UNLC25	Weight: 2.18 kg Capacity: 25 litres	External Size: 500 x 300 x 274mm Internal Size: 422 x 227 x 239mm
UNLC45	Weight: 3.26 kg Capacity: 45 litres	External Size: 600 x 400 x 275mm Internal Size: 528 x 328 x 238mm
UNLC60	Weight: 4.20 kg Capacity: 60 litres	External Size: 600 x 400 x 374mm Internal Size: 527 x 315 x 335mm

SUPPLEMENTARY PRODUCTS

CODE	DESCRIPTION	QUANTITY
UNLC-Y	Sealing pin for packing group ii and iii	100
UNLC-X	Sealing strip for packing group i, ii and iii	100

Transportation Accessories

A range of accessories to guarantee the identification and security of transported medical instruments and equipment.

Transportation box labels

LF

A selection of small and medium-sized printed labels designed to be inserted within the plastic pockets of BTB and RTB transportation boxes. Our biohazard labels conform to the requirements for a 6.2 hazard (infectious substances) as stated within ANNEX A (sub-section 5.2.2.2.2) of the European Agreement concerning the International Carriage of Dangerous Goods by Road. Laminated labels with the words 'USED MEDICAL DEVICE' are also available for the transportation of used medical devices en route for decontamination.



CLASS 6.2 BIOHAZARD LABEL

CODE	DESCRIPTION	SIZE	QUANTITY
BTBLABEL001	Non-laminated label	125 x 40mm	20
BTBLABEL002	Non-laminated label	150 x 70mm	20
BTBLABEL003	Laminated label	125 x 40mm	20
BTBLABEL004	Laminated label	150 x 70mm	20



USED MEDICAL DEVICE LABEL

CODE	DESCRIPTION	SIZE	QUANTITY
BTB232SML	Laminated label	125 x 40mm	100
BTB232LRG	Laminated label	150 x 70mm	100

Plastic cable ties



COLOUR CABLE TIES

CODE	DESCRIPTION	SIZE	QUANTITY
TIE001	Cable tie blue	144mm x 3.6mm	100
TIE007	Cable tie large yellow	200mm x 4.8mm	100
TIE008	Cable tie large green	300mm x 4.8mm	100
TIE009	Re-sealable cable tie with dual closure design. Not to be used in an autoclave.	120mm x 3.5mm	1000



NUMBERED CABLE TIES

CODE	DESCRIPTION	SIZE	QUANTITY
TIE003	Individually numbered tie (green)	230mm x 18mm	100
TIE005	Individually numbered tie (red)	230mm x 18mm	100



SECURITY CABLE TIES

CODE	DESCRIPTION	SIZE	QUANTITY
TIE010	Red security tag tie. Not to be used in an autoclave.	300mm x 4.8mm	100
BTBREDTIE	Red security tie. This is also available in yellow (BTBYELTIE) and green (BTBGRNTIE)	145mm x 3.6mm	100



CLEAR CABLE TIES

CODE	DESCRIPTION	SIZE	QUANTITY
TIE002	Cable tie clear	144mm x 2.5mm	100
TIE004	Cable tie clear (with additional thickness)	300mm x 3.6mm	100
TIE006	Secur-a-tie clear	125mm x 1.5mm	500

BTB Trolley



W LF

CODE	DESCRIPTION
BTB004TROL	Trolley to fit BTB004/RTB004
BTB005TROL	Trolley to fit BTB005/RTB005
BTB006TROL	Trolley to fit BTB006/RTB006

BTB Opener



W LF

CODE	DESCRIPTION
BTBOPENER	Stainless steel opener suitable for all BTB and RTB boxes

"... Clinipak has really supported UCLH with brilliant customer services and deliveries always on time."

Joe Curran
Materials Manager
University College London Hospital NHS Foundation Trust



Transportation pouches

An innovative solution to assist when the immediate cleaning of surgical instruments is not possible.

Humipak



Humipak is a patented method to prevent any protein from drying out on instruments before subsequent cleaning.

Each Humipak pouch comprises of two layers of film with a highly absorbent material inside its interior to ensure moisture is retained once sealed. This product is implemented most effectively when the specified amount of water is added into the pouch until the instruments inside are manually cleaned in the CSSD. The range consists of three large pouches to hold baskets of the following sizes: half DIN, DIN and double DIN. A small edition of the pouch has been produced to accommodate smaller baskets. The pouches demonstrate the capacity to retain a moist atmosphere for a seven day period.



CODE	DESCRIPTION	QUANTITY	SIZE
HPSS2035	Small	100/carton	200 x 350mm
HPSS4243	Half 'DIN'	100/carton	420 x 430mm
HPSS4267	Full 'DIN'	100/carton	420 x 670mm
HPSS6577	Double 'DIN'	50/carton	650 x 770mm



Storage

pages 49-52



Storage Systems

A selection of products to support SSD units with the storage of sterile sets. This has been achieved by means of effective storage systems and transport trolleys.

Belintra UBeFlex storage system



The UBeFlex® modular storage system was developed with very quick assembly, functionality, durability and stability in mind. The result is a unique storage system made from anodised aluminium and other high-quality, anti-static synthetic materials, which means it can be used in both sterile and non-sterile environments.

Thanks to the modular construction, in D-type or W-type columns, which may or may not be combined with shelves, it can be used in various departments of the healthcare institution, from CSA through pharmacy to kitchen and cold storage. Because of the unique design, it is an ergonomic and hygienic storage solution and can also be used as a "pass-through" rack.

Please contact a member of our sales team for further product information.



Belintra Stainless Steel Trolley



Belintra stainless steel trolleys are resistant to washing tunnels, fitted with a lock with breakable seal, and produced in accordance with Euronorm standard 600 x 400 mm or DIN standard 600 x 300 mm. This Belintra transport trolley has been designed to transport stored sterile sets between external sterilisation units and hospitals (and also from hospitals to central decontamination units). These trolleys are also compatible with cart washers.



- 270°C swivel doors
- Removeable runners
- 125mm diameter wheels
- 160mm diameter castors

Belintra UBeFlex KD Transport Trolley



The Belintra Uflex KD Transport trolley is produced with anodised aluminium, equipped with nylon 'no-noise' runners and designed for the internal transport of sterile and non-sterile goods. A polypropylene cover can be introduced to facilitate the use of the trolley in sterile environments.

The trolley is easy to adapt, demonstrates highly durability and exhibits non-marking swivel wheels which are easily maneuvered and can be locked into place.



Transportation trays and bags

A range of products to ensure that all sterile sets can be further protected from dust and rigorous handling.

Belintra transportation trays



Stainless steel transport trays provide protection to packs during transportation.



CODE	SIZE
BELTRY010	600 x 400 x 100mm
BELTRY013	500 x 540 x 70mm

Clinipak transportation trays

LF

Clinipak transportation trays reduce the damage that wrapped packs can endure as a result of poor handling and shelving during the transport and storage stages of the decontamination process. This product is available in sealed-edge polypropylene (correx) or cardboard. Both variants can be used alone or with a strap and two transportation trays can be interlocked to form a lidded box.

The unique sealed polypropylene board in particular enables the trays to be processed through a decontamination wash cycle for repeated use.



POLYPROPYLENE TRANSPORTATION TRAYS

CODE	SIZE
TTRAY002	500 x 270 x 70mm

CARDBOARD TRANSPORTATION TRAYS

CODE	SIZE
TTRAY001	500 x 270 x 70mm
TTRAY003	500 x 540 x 70mm



CLINIPAK TRAY STRAPS

CODE	SIZE	QUANTITY
STR002	65cm / Yellow	10
STR004	140cm / Blue	10

Transportation bags

2 LF

Clinipak transportation bags are manufactured with high-quality polyethylene to protect sterile packs from rigorous handling and provide an effective dust cover. All bags exhibit tamper-proof seals, tears which can be simply opened, and 'half moon' ventilation holes to prevent condensate from spreading in a humid environment.



CODE	SIZE	QUANTITY
TBAG001	380 x 760mm	250
TBAG002	585 x 765mm	200
TBAG003	380 x 710mm	250
TBAG004	385 x 490mm	250



Use

pages 53-56

Theatre sundries

A selection of sundry products for use in Theatres ranging from hand hygiene products to apparel and accessories.

Medica scrub brushes

A reputable selection of three disposable pre-operative scrub brushes which are available sterile or pre-soaped. All Medica scrub brushes are suitable for use in Operating Theatres and other clinical areas for cleaning hands and skin.

CODE	DESCRIPTION	QUANTITY	PROPERTIES
MED001	Dry/sterile	15 boxes of 25 brushes	2 CE LF
MED002	Chlorhexidine	15 boxes of 25 brushes	2 LF
MEDDRYSML25 (previously FSL061)	Dry/Sterile	1 pack of 25 brushes	2 CE LF
MEDCHLSML25 (previously FSL070)	Chlorhexidine	1 pack of 25 brushes	2 LF



O.R. soaker sheets



Clinipak O.R. soakers are designed to be used in Operating Theatres to absorb fluid spillages during procedures. The non-slip backing prevents movement when laid on the floor and the upper surface absorbs between 1 and 3 litres per m². This product is of particular use in Urology and Obstetrics cases.



CODE	DESCRIPTION	SIZE
ORSOAK1000	Absorption capacity of 1 litre/m ²	60 x 100cm
ORSOAK3000	Absorption capacity of 3 litre/m ²	60 x 100cm

Foam cubes and sheets

S ~~2~~ CE LF Custom sizes available

Foam cubes and sheets are produced with white standard foam and supplied in any size and unsterile. Our flexible production methods ensure that custom sizes can be provided on request.



CODE	SIZE	QUANTITY
FOB009	50 x 38 x 38mm	1000/ctn
FOB010	50 x 38 x 38mm	1000/ctn
FOB012	38 x 38 x 38mm	1000/ctn
FOB046	50 x 50 x 50mm	1000/ctn
FOS050	254 x 254 x 25mm	1000/ctn

Gown ties

S ~~2~~ LF

An effective solution to hold gown ties together and facilitate aseptic donning procedures. The unique pattern allows for ties to be inserted into both ends or through its centre for simple tie shortening. This product is of particular use in Laundry Departments to ensure that all re-usable gown ties are well organised.



CODE	DESCRIPTION	QUANTITY
GOW001	Combination gown ties	7000
GOW002	Crocodile cut for micro fibre gowns	7000

Latex-free silicone bands

S EO ~~2~~ LF

Clinipak latex-free silicone bands are increasingly popular due to higher awareness of latex allergies. These bands can be used in various situations in which natural latex can cause problems (e.g. hand surgery, securing instruments, flap retention, bundling of arterial forceps etc.)



CODE	SIZE	QUANTITY
BAND001	60 x 1.5 x 3mm	1000
BAND002	75 x 1.5 x 3mm	800

Clogs



MedisMore theatre clogs have been manufactured with leading edge technology and designed for maximum comfort. This entire selection of clogs is supplied without upper holes for ultimate protection. Available sizes range from UK size 3.5 (EUR size 36) to UK size 11 (EUR size 46).

Heel straps provided on request – add 'S' to the product code.



CODE	DESCRIPTION	SIZE	QUANTITY
CLOGNH(36-46)DG	Dark green	3.5 - 11 / 36 - 46	1 pair
CLOGNH(36-46)LB	Light blue	3.5 - 11 / 36 - 46	1 pair
CLOGNH(36-45)PK	Pink	3.5 - 10 / 36 - 45	1 pair
CLOGNH(36-45)DB	Dark blue	3.5 - 10 / 36 - 45	1 pair
CLOGNH(36-45)LL	Lilac	3.5 - 10 / 36 - 45	1 pair

• MedisMore Theatre Clogs are compliant with the requirements of the European Directive 89/686/EEC relating to personal protective equipment and are CE-marked on the sole of the clog.

Clog cart

CODE	DESCRIPTION	SIZE
CLOGCART	A movable trolley which can hold a maximum of 30 pairs of clogs. This clog cart maximises space and improves health and safety as clogs are well organised. Flat packed for easy self-assembly.	90 x 60cm



LF Clog rack

CODE	DESCRIPTION	SIZE
CLOGRACK	A racking system that facilitates the arrangement of 15 pairs of clogs. This clog rack is especially useful in changing or operating rooms.	90 x 60cm





Please contact a member of our sales team if you are looking for trusted solutions to support your decontamination process.



Email: sales@clinipak.co.uk
Telephone: +44 (0) 1628 810626
Fax: +44 (0) 1628 810760

Index (A-Z)

Absorbent and protective mats	33	Gown ties	55
Aluminum trays	15	Healthmark trays	16
Belintra premium baskets	12	Humipak	48
Belintra stainless steel trolley	50	Identification tags	21
Belintra Steri-ID range	12	Instrument locks	18
Belintra transportation trays	51	Latex-free tapes	28
Belintra UFlex KD transport trolley	51	Latex-free silicone bands	55
Belintra UFlex shelving system	50	Lidded crates (UN3291)	45
Blue brushes	08	Marker pens	21
Bowie and Dick packs	40	Marking tape	20
BTB Opener	47	Medica scrub brushes	54
BTB Trolley	47	Microclene	06
CBM baskets	13	Microclene dispenser	07
CBM container filters	42	Millenium liners	32
CBM containers	41	Nic-nac baskets	11
Chessboard instrument protector	30	Numbered instrument locks	22
Choice bags	26	O.R. soaker sheets	54
Choice dry bags	26	Paper multipouches	31
Choice wrap	25	Plastic cable ties	46
Choice wrap dry	25	Plastic scope trays	14
Chux low-linting absorbant towels (now known as Veraclean)	18	Plastic sterilisation trays	16
Clinipak 340 wrap	24	Point protectors	30
Clinipak corners	29	polySteribox®	42
Clinipak premium baskets	10	Security ties for tray tags	20
Clinipak standard pouches & reels	34	Self-adhesive alert labels	22
Clinipak transportation trays	52	Silicone inserts	14
Clog cart	56	Silicone matting	31
Clog rack	56	Split rings	20
Clogs	56	Stainless steel scope baskets	14
Container padlocks	40	Standard autoclave tapes	28
Crepe and plain paper wrap	24	Soaker sheets	33
Cut corner liners	33	Steriking premium pouches & reels	36
Delicate wash tray baskets	15	Sterilisation paper bags	38
Diathermy cleaning brushes	08	Tape dispensers	28
Dispensers & cutters	37	Transparent point protectors	30
Dry transportation cover wraps	27	Transportation cover bags	26
Endoscopy brushes	08	Transportation cover wrap	27
Flat base economy baskets	10	Transportation bags	52
Flexistem	06	Transportation boxes (UN3291)	44
Flexistem dispenser	07	Transportation box labels	45
Foam instrument pouches	31	Tray anchors	29
Foam cubes and sheets	55	Tray dividers	11
		Tray liners	32
		Tray tags	19
		Trinket baskets	11
		Tyvek® pouches & reels	35
		Washer tags	19
		Wavy flexistem	06
		WIPAK multi-seal test kit	41
		WI SCAN	41
		Wire mesh economy baskets	10

Clinipak opening hours:
Monday to Friday (8:00am – 5:00pm)

Please visit our website to find additional
product information and special offers.

www.clinipak.co.uk

CLINIPAK

Product Catalogue 2015
Issue 4.6



Tel. +44 (0) 1628 810626
Fax. +44 (0) 1628 810760
sales@clinipak.co.uk
www.clinipak.co.uk

Clinipak Limited
Beech House
Knaves Beech Business Centre
Loudwater
High Wycombe
Buckinghamshire
HP10 9SD
United Kingdom

