

KEY SURGICAL INSTRUMENT CLEANING

The Key Surgical Instrument Cleaning line is designed to provide you with the brushes needed to clean your surgical instruments. Our wide variety of cleaning brushes ensures you and your team have the right brush at the right time.

With bristles available in nylon, brass, polypropylene or stainless steel, we also provide an incredible variety of lengths and diameters to fit your needs.

For organization of your cleaning brushes, our brush organizer is perfect for hanging in the department. Our stainless steel organizer is easy to keep clean and helps your team choose the right brush in a more efficient manner.

With our vast assortment of brushes - from channel cleaning to scope cleaning to sterilizer cleaning - we can help you match the right brush to the right cleaning task at hand.

Channel Cleaning Brushes

Pipe Cleaners

Toothbrush-Style Cleaning Brushes

General Instrument Cleaning Brushes

Scope Cleaning Brushes

Brush Assortments


Sterilizer Cleaning Brushes

Brush Organizer

French Catheter Scale

PLEASE NOTE: Unless specified, all products are not made with natural rubber latex.





BRUSH TIP DETAILS SHOWN ACTUAL SIZE


Instrument Cleaning – Channel Cleaning Brushes

CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. FR indicates equivalent French catheter size.

	Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A	BR-06-039	.039"	1.00mm	3	.25"	6"	2
	BR-06-039-50						50
B	BR-08-039	.039"	1.00mm	3	.50"	8"	2
	BR-08-039-50						50
	BR-18-039	.039"	1.00mm	3	.25"	18"	2
	BR-18-039-50						50
C	BR-06-062	.062"	1.58mm	5	.50"	6"	2
	BR-06-062-50						50
	BR-18-062	.062"	1.58mm	5	.50"	18"	2
	BR-18-062-50						50
D	BR-08-070	.070"	1.78mm	5	.50"	8"	2
	BR-08-070-50						50
	BR-18-070	.070"	1.78mm	5	.50"	18"	2
	BR-18-070-50						50
	BR-24-070	.070"	1.78mm	5	.50"	24"	2
	BR-24-070-50						50
E	BR-06-079	.079"	2.00mm	6	.50"	6"	2
	BR-06-079-50						50
F	BR-08-079	.079"	2.00mm	6	.75"	8"	2
	BR-08-079-50						50
	BR-12-079	.079"	2.00mm	6	.75"	12"	2
	BR-12-079-50						50
G	BR-24-080	.080"	2.03mm	6	1"	24"	2
	BR-24-080-50						50
H	BR-08-087	.087"	2.20mm	7	.50"	8"	2
	BR-08-087-50						50
	BR-18-087	.087"	2.20mm	7	.50"	18"	2
	BR-18-087-50						50
I	BR-08-092	.092"	2.34mm	7	.50"	8"	2
	BR-08-092-50						50
J	BR-24-098	.098"	2.50mm	8	.50"	24"	2
	BR-24-098-50						50
K	BR-08-105	.105"	2.67mm	8	1"	8"	2
	BR-08-105-50						50
	BR-18-105	.105"	2.67mm	8	1"	18"	2
	BR-18-105-50						50
L	BR-08-118	.118"	3.00mm	9	1.50"	8"	2
	BR-08-118-50						50
M	BR-12-118	.118"	3.00mm	9	2.75"	12"	2
	BR-12-118-50						50





BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes



CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. FR indicates equivalent French catheter size.

	Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A	BR-18-125	.125"	3.18mm	10	2"	18"	2
	BR-18-125-50						50
B	BR-08-131	.131"	3.33mm	10	.75"	8"	2
	BR-08-131-50						50
C	BR-08-144	.144"	3.67mm	11	.75"	8"	2
	BR-08-144-50						50
D	BR-08-158	.158"	4.00mm	12	.75"	8"	2
	BR-08-158-50						50
	BR-12-158	.158"	4.00mm	12	.75"	12"	2
	BR-12-158-50						50
	BR-18-158	.158"	4.00mm	12	.75"	18"	2
	BR-18-158-50						50
E	BR-12-197	.197"	5.00mm	15	3"	12"	2
	BR-12-197-50						50
	BR-16-197	.197"	5.00mm	15	3"	16"	2
	BR-16-197-50						50
	BR-18-197	.197"	5.00mm	15	3"	18"	2
	BR-18-197-50						50
F	BR-24-197	.197"	5.00mm	15	3"	24"	2
	BR-24-197-50						50
	BR-24-200	.200"	5.08mm	15	2"	24"	2
	BR-24-200-50						50
G	BR-12-236	.236"	6.00mm	18	2"	12"	2
	BR-12-236-50						50
	BR-16-236	.236"	6.00mm	18	2"	16"	2
	BR-16-236-50						50
	BR-24-236	.236"	6.00mm	18	2"	24"	2
	BR-24-236-50						50



A



B



C



D



E



F



G

BRUSH TIP DETAILS SHOWN ACTUAL SIZE


Instrument Cleaning – Channel Cleaning Brushes



CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. FR indicates equivalent French catheter size.

	Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A	BR-06-250	.250"	6.35mm	19	2"	6"	2
	BR-06-250-50						50
	BR-08-250	.250"	6.35mm	19	2"	8"	2
	BR-08-250-50						50
	BR-18-250	.250"	6.35mm	19	2"	18"	2
	BR-18-250-50						50
B	BR-16-276	.276"	7.00mm	21	3"	16"	2
	BR-16-276-50						50
	BR-24-276	.276"	7.00mm	21	3"	24"	2
	BR-24-276-50						50
C	BR-24-312	.312"	8.00mm	24	3"	24"	2
	BR-24-312-50						50
D	BR-06-375	.375"	9.52mm	29	1.5"	6"	2
	BR-06-375-50						50
E	BR-24-375	.375"	9.52mm	29	3"	24"	2
	BR-24-375-50						50
F	BR-16-393	.393"	10.00mm	30	3"	16"	2
	BR-16-393-50						50
	BR-18-393	.393"	10.00mm	30	3"	18"	2
	BR-18-393-50						50
	BR-24-393	.393"	10.00mm	30	3"	24"	2
	BR-24-393-50						50
G	BR-24-425	.425"	10.80mm	32	3"	24"	2
	BR-24-425-50						50



BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes



CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. FR indicates equivalent French catheter size.

	Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A	BR-08-440	.440"	11.18mm	34	3"	8"	2
	BR-08-440-50						50
	BR-24-440	.440"	11.18mm	34	3"	24"	2
	BR-24-440-50						50
B	BR-16-443	.443"	11.25mm	34	3"	16"	2
	BR-16-443-50						50
	BR-24-443	.443"	11.25mm	34	3"	24"	2
	BR-24-443-50						50
C	BR-16-472	.472"	12.00mm	N/A	3"	16"	2
	BR-16-472-50						50
D	BR-06-500	.500"	12.70mm	N/A	2"	6"	2
	BR-06-500-50						50
E	BR-16-591	.591"	15.00mm	N/A	2"	16"	2
	BR-16-591-50						50
	BR-18-591	.591"	15.00mm	N/A	2"	18"	2
	BR-18-591-50						50

A LITTLE BIT ABOUT...

Brian O'Connell

Job Title: **President and Chief Operating Officer**

Who inspires you? **Great scientific thinkers; like Einstein and Stephen Hawking.** What is your favorite thing to watch on TV? **Anything on the Science Channel, Smithsonian Channel, or any episode of The Sopranos.** If you were to perform in the circus, what would you perform? **Some sort of aerial act for sure.** If you could be a teacher of something, what would it be? **Astronomy. I find the science and physics of the universe incredibly interesting. Quantum theory and black holes – remarkable.** Who was your favorite band growing up? And your favorite band now? **Growing up – The Who. Now, pretty much anything AOC (adult oriented contemporary). Groups like Muse, Bon Iver, Daft Punk as well as Roxy Music, David Bowie, and Split Enz.** What do you do to relax? **Play golf.**



BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes




CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end.

Number	Diameter inches	Diameter mm	Bristle Area Length	Length	Quantity/pkg
A BR-08-630	.630"	16.00mm	2"	8"	2
BR-08-630-50					50
B BR-16-787	.787"	20.00mm	2"	16"	2
BR-16-787-50					50
C BR-18-880	.880"	22.35mm	3.50"	18"	2
BR-18-880-50					50
D BR-24-880	.880"	22.35mm	4.25"	24"	2
BR-24-880-50					50

FUN FACT:

The three women who started Key Surgical in 1988 had years of experience in both the nursing profession and the medical supply industry. They were motivated by the growing need for a single-source supplier that would provide excellent products and equally excellent customer service. One of the founders is Gail Milchman, the mother of our CEO, Scot Milchman.



Instrument Cleaning – Channel Cleaning Brushes



CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end (pull through style excluded). FR indicates equivalent French catheter size.

	Number	Details	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A	BR-3044	nylon	.020"	.50mm	2	.25"	3"	10
	BR-3045	nylon	.020"	.50mm	2	.25"	8"	10
B	BR-3047	nylon	.035"	.90mm	3	.5"	4"	10
C	BR-4120	nylon, tapered	.250" to .600"	6.35mm to 15.24mm	N/A	2"	18"	2
	BR-4120-50							50
D	BR-4140	nylon, centered, pull through	.500"	12.70mm	N/A	3"	24"	6



BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes



CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. FR indicates equivalent French catheter size.

	Number	Description	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A	BR-8110	curved bristle tip	.394"	10mm	29	1.75"	4.6"	3
	BR-8110-50							50
B	BR-8115	curved bristle tip	.551"	13.99mm	NA	3.625"	8"	3
	BR-8115-50							50
C	BR-8210	tapered fan bristle tip	.187" to .500"	4.75mm to 12.70mm	NA	3"	7"	3
	BR-8210-50							50
D	BR-8310	straight bristle tip	.500"	12.70mm	NA	2.5"	12"	3
	BR-8310-50							50
E	BR-8315	straight bristle tip	.500"	12.70mm	NA	2.5"	12"	3
	BR-8315-50							50
F	BR-8410	fan bristle tip	.500"	12.70mm	NA	3"	6"	3
	BR-8410-50							50

DID YOU KNOW?

The Key Surgical Community Giving Committee seeks out volunteer opportunities for all Key Surgical employees. We partner with local organizations that provide services and assistance to low-income families, people with disabilities, children, the homeless, and those suffering from terminal illness.



Instrument Cleaning – Channel Cleaning Brushes




BONE REAMER BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a rigid stainless steel handle with loop end. FR indicates equivalent French catheter size.

	Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A	BR-31-118	.118"	3.00mm	9	3"	31"	3
B	BR-31-197	.197"	5.00mm	15	3"	31"	3
C	BR-31-276	.276"	7.00mm	21	3"	31"	3





BRUSH TIP DETAILS SHOWN ACTUAL SIZE


Instrument Cleaning – Channel Cleaning Brushes



LUMEN GUARD™ CLEANING BRUSHES

Channel cleaning brushes featuring rigid, antimicrobial nylon bristles and a protective tip on a twisted stainless steel handle with loop end. Inseparable, protective tip created using patented Z-tip technology. FR indicates equivalent French catheter size.

	Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A	BZ-10-050 BZ-10-050-50	.050"	1.27mm	4	1"	10"	2 50
B	BZ-12-062 BZ-12-062-50	.062"	1.58mm	5	1"	12"	2 50
C	BZ-18-062 BZ-18-062-50	.062"	1.58mm	5	2"	18"	2 50
D	BZ-18-079 BZ-18-079-50	.079"	2.00mm	6	2"	18"	2 50
E	BZ-12-087 BZ-12-087-50	.087"	2.20mm	7	1"	12"	2 50
F	BZ-12-098 BZ-12-098-50	.098"	2.50mm	8	1"	12"	2 50
G	BZ-18-105 BZ-18-105-50	.105"	2.67mm	8	2"	18"	2 50
H	BZ-08-125 BZ-08-125-50	.125"	3.18mm	10	.75"	8"	2 50
I	BZ-12-125 BZ-12-125-50	.125"	3.18mm	10	1"	12"	2 50
J	BZ-18-125 BZ-18-125-50	.125"	3.18mm	10	2"	18"	2 50
	BZ-24-125 BZ-24-125-50	.125"	3.18mm	10	2"	24"	2 50
K	BZ-18-131 BZ-18-131-50	.131"	3.33mm	10	2"	18"	2 50
L	BZ-18-144 BZ-18-144-50	.144"	3.37mm	11	2"	10"	2 50
M	BZ-12-158 BZ-12-158-50	.158"	4.00mm	12	1"	12"	2 50
N	BZ-18-158 BZ-18-158-50	.158"	4.00mm	12	2"	18"	2 50



BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes



LUMEN GUARD™ CLEANING BRUSHES (CONTINUED)

Channel cleaning brushes featuring rigid, antimicrobial nylon bristles and a protective tip on a twisted stainless steel handle with loop end. Inseparable, protective tip created using patented Z-tip technology. FR indicates equivalent French catheter size.

	Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A	BZ-18-197 BZ-18-197-50	.197"	5.00mm	15	2"	18"	2 50
B	BZ-18-236 BZ-18-236-50	.236"	6.00mm	18	2"	18"	2 50
C	BZ-16-250 BZ-16-250-50 BZ-18-250 BZ-18-250-50	.250" .250"	6.35mm 6.35mm	19 19	2" 2"	16" 18"	2 50 2 50
D	BZ-18-312 BZ-18-312-50	.312"	8.00mm	24	2"	18"	2 50
E	BZ-18-354 BZ-18-354-50	.354"	9.00mm	27	2"	18"	2 50
F	BZ-18-393 BZ-18-393-50	.393"	10.00mm	30	2"	18"	2 50
G	BZ-18-433 BZ-18-433-50	.433"	11.00mm	33	2"	18"	2 50
H	BZ-18-472 BZ-18-472-50	.472"	12.00mm	N/A	2"	18"	2 50
I	BZ-18-1000 BZ-18-1000-50	1"	25.40mm	N/A	2"	18"	2 50
J	BZ-18-1543 BZ-18-1543-50	.157" to .433"	4mm to 11mm	N/A	2"	18"	2 50
K	BZ-06-2550 BZ-06-2550-50	.250" to .500"	6.35mm to 12.7mm	N/A	.75"	6"	2 50



A



B



C



D



E



F



G



H



I

BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes



FAN TIP BRUSHES

Channel cleaning brushes featuring a fan tip to assist with cleaning closed-end channels. Antimicrobial nylon bristles on a twisted stainless steel handle with loop end. FR indicates equivalent French catheter size.

	Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A	FT-12-197	.197"	5.00mm	15	2"	12"	2
	FT-12-197-50						50
	FT-18-197	.197"	5.00mm	15	2"	18"	2
	FT-18-197-50						50
	FT-24-197	.197"	5.00mm	15	3"	24"	2
	FT-24-197-50						50
B	FT-12-250	.250"	6.35mm	19	2"	12"	2
	FT-12-250-50						50
	FT-18-250	.250"	6.35mm	19	3"	18"	2
	FT-18-250-50						50
	FT-24-250	.250"	6.35mm	19	3"	24"	2
	FT-24-250-50						50
C	FT-12-276	.276"	7.00mm	21	2"	12"	2
	FT-12-276-50						50
	FT-18-276	.276"	7.00mm	21	3"	18"	2
	FT-18-276-50						50
	FT-24-276	.276"	7.00mm	21	3"	24"	2
	FT-24-276-50						50
D	FT-12-375	.375"	9.52mm	29	2"	12"	2
	FT-12-375-50						50
E	FT-18-375	.375"	9.52mm	29	3"	18"	2
	FT-18-375-50						50
	FT-24-375	.375"	9.52mm	29	3"	24"	2
	FT-24-375-50						50
F	FT-12-440	.440"	11.18mm	34	2"	12"	2
	FT-12-440-50						50
G	FT-18-440	.440"	11.18mm	34	3"	18"	2
	FT-18-440-50						50
	FT-24-440	.440"	11.18mm	34	3"	24"	2
	FT-24-440-50						50
H	FT-12-500	.500"	12.70mm	N/A	2"	12"	2
	FT-12-500-50						50
I	FT-18-500	.500"	12.70mm	N/A	3"	18"	2
	FT-18-500-50						50
	FT-24-500	.500"	12.70mm	N/A	3"	24"	2
	FT-24-500-50						50



BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes

ACRYLIC TIP BRUSHES

Channel cleaning brushes featuring a protective acrylic tip for an extra level of protection during cleaning. Antimicrobial nylon bristles on a flexible handle. No loop end, may be used as a pull through brush. FR indicates French catheter size.

	Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A	AT-12-062	.062"	1.58mm	5	1"	12"	2
	AT-12-062-50						50
	AT-18-062	.062"	1.58mm	5	1"	18"	2
	AT-18-062-50						50
	AT-24-062	.062"	1.58mm	5	1"	24"	2
	AT-24-062-50						50
B	AT-12-087	.087"	2.20mm	7	1"	12"	2
	AT-12-087-50						50
	AT-18-087	.087"	2.20mm	7	1"	18"	2
	AT-18-087-50						50
	AT-24-087	.087"	2.20mm	7	1"	24"	2
	AT-24-087-50						50
C	AT-12-105	.105"	2.67mm	8	1"	12"	2
	AT-12-105-50						50
	AT-18-105	.105"	2.67mm	8	1"	18"	2
	AT-18-105-50						50
	AT-24-105	.105"	2.67mm	8	1"	24"	2
	AT-24-105-50						50
D	AT-12-118	.118"	3.00mm	9	1"	12"	2
	AT-12-118-50						50
	AT-18-118	.118"	3.00mm	9	1"	18"	2
	AT-18-118-50						50
	AT-24-118	.118"	3.00mm	9	1"	24"	2
	AT-24-118-50						50
E	AT-12-125	.125"	3.18mm	10	1"	12"	2
	AT-12-125-50						50
	AT-18-125	.125"	3.18mm	10	1"	18"	2
	AT-18-125-50						50
	AT-24-125	.125"	3.18mm	10	1"	24"	2
	AT-24-125-50						50
F	AT-12-131	.131"	3.33mm	10	1"	12"	2
	AT-12-131-50						50
	AT-18-131	.131"	3.33mm	10	1"	18"	2
	AT-18-131-50						50
	AT-24-131	.131"	3.33mm	10	1"	24"	2
	AT-24-131-50						50





A



B



C



D



E



F

BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes

ACRYLIC TIP BRUSHES (CONTINUED)

Channel cleaning brushes featuring a protective acrylic tip for an extra level of protection during cleaning. Antimicrobial nylon bristles on a flexible handle. No loop end, may be used as a pull through brush. FR indicates French catheter size.

	Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A	AT-12-158	.158"	4.00mm	12	1.5"	12"	2
	AT-12-158-50						50
	AT-18-158	.158"	4.00mm	12	1.5"	18"	2
	AT-18-158-50						50
	AT-24-158	.158"	4.00mm	12	1.5"	24"	2
	AT-24-158-50						50
B	AT-12-197	.197"	5.00mm	15	2"	12"	2
	AT-12-197-50						50
	AT-18-197	.197"	5.00mm	15	2"	18"	2
	AT-18-197-50						50
	AT-24-197	.197"	5.00mm	15	2"	24"	2
	AT-24-197-50						50
C	AT-12-250	.250"	6.35mm	19	2"	12"	2
	AT-12-250-50						50
	AT-18-250	.250"	6.35mm	19	2"	18"	2
	AT-18-250-50						50
	AT-24-250	.250"	6.35mm	19	2"	24"	2
	AT-24-250-50						50
D	AT-12-276	.276"	7.00mm	21	2"	12"	2
	AT-12-276-50						50
	AT-18-276	.276"	7.00mm	21	2"	18"	2
	AT-18-276-50						50
	AT-24-276	.276"	7.00mm	21	2"	24"	2
	AT-24-276-50						50
E	AT-12-375	.375"	9.52mm	29	3"	12"	2
	AT-12-375-50						50
	AT-18-375	.375"	9.52mm	29	3"	18"	2
	AT-18-375-50						50
	AT-24-375	.375"	9.52mm	29	3"	24"	2
	AT-24-375-50						50
F	AT-12-440	.440"	11.18mm	34	3"	12"	2
	AT-12-440-50						50
	AT-18-440	.440"	11.18mm	34	3"	18"	2
	AT-18-440-50						50
	AT-24-440	.440"	11.18mm	34	3"	24"	2
	AT-24-440-50						50





A



B



C



D



E



F



G

BRUSH TIP DETAILS SHOWN ACTUAL SIZE


Instrument Cleaning – Channel Cleaning Brushes



POLYPROPYLENE BRUSHES

Channel cleaning brushes featuring rigid polypropylene bristles on a twisted stainless steel handle with loop end. FR indicates equivalent French catheter size.

	Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A	PP-16-039	.039"	1.00mm	3	1"	16"	2
	PP-16-039-50						50
	PP-22-039	.039"	1.00mm	3	1"	22"	2
	PP-22-039-50						50
B	PP-16-062	.062"	1.58mm	5	1"	16"	2
	PP-16-062-50						50
	PP-22-062	.062"	1.58mm	5	1"	22"	2
	PP-22-062-50						50
C	PP-16-087	.087"	2.20mm	7	1"	16"	2
	PP-16-087-50						50
	PP-22-087	.087"	2.20mm	7	1"	22"	2
	PP-22-087-50						50
D	PP-16-105	.105"	2.67mm	8	2"	16"	2
	PP-16-105-50						50
	PP-22-105	.105"	2.67mm	8	2"	22"	2
	PP-22-105-50						50
E	PP-16-125	.125"	3.18mm	10	2"	16"	2
	PP-16-125-50						50
	PP-22-125	.125"	3.18mm	10	2"	22"	2
	PP-22-125-50						50
F	PP-16-158	.158"	4.00mm	12	2"	16"	2
	PP-16-158-50						50
	PP-22-158	.158"	4.00mm	12	2"	22"	2
	PP-22-158-50						50
G	PP-18-197	.197"	5.00mm	15	2"	18"	2
	PP-18-197-50						50
	PP-22-197	.197"	5.00mm	15	2"	22"	2
	PP-22-197-50						50



H



I



J



K



L



M

(PULL THROUGH)

BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes



POLYPROPYLENE BRUSHES (CONTINUED)

Channel cleaning brushes featuring rigid polypropylene bristles on a twisted stainless steel handle with loop end (pull through styles excluded). FR indicates equivalent French catheter size.

	Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
H	PP-18-250	.250"	6.35mm	19	2"	18"	2
	PP-18-250-50						50
	PP-24-250	.250"	6.35mm	19	2"	24"	2
	PP-24-250-50						50
I	PP-18-393	.393"	10.00mm	30	2"	18"	2
	PP-18-393-50						50
J	PP-18-591	.591"	15.00mm	N/A	2"	18"	2
	PP-18-591-50						50
K	PP-18-750	.750"	19.05mm	N/A	6"	18"	2
	PP-18-750-50						50
L	PT-28-118	.118"	3.00mm	9	2"	28"	2
	PT-28-118-50 (pull through)						50
M	PT-28-197	.197"	5.00mm	15	2"	28"	2
	PT-28-197-50 (pull through)						50



DIAMETERS SHOWN ACTUAL SIZE

Instrument Cleaning – Pipe Cleaners



PIPE CLEANERS - 50' ROLLS

Available in either cotton or polyester. Rolls allow you to cut to the size you need, minimizing waste. Single use.

	Number	Material	Diameter inches	Diameter mm	Quantity/pkg
A	PC-100C	cotton	.125"	3.18mm	1 roll
B	PC-200C	cotton	.250"	6.35mm	1 roll
C	PC-300P	polyester	.125"	3.18mm	1 roll
D	PC-400P	polyester	.250"	6.35mm	1 roll

PIPE CLEANERS - 12" LENGTHS

E	PC-112C	cotton	.125"	3.18mm	100/pkg
F	PC-212C	cotton	.250"	6.35mm	100/pkg
G	PC-312P	polyester	.125"	3.18mm	100/pkg
H	PC-412P	polyester	.250"	6.35mm	100/pkg

DOUBLE CLEANER

Soft cotton and nylon cleaner. Single use.

	Number	Diameter inches	Diameter mm	Length	Quantity/pkg
I	PC-503	.118"	3mm	6"	100/pkg

FLEECE CLEANER

Fleece cleaner for extra delicate items. Single use.

J	PC-512	.500"	12.7mm	12"	50/pkg
----------	--------	-------	--------	-----	--------

FLEX CLEANER

Polyester cleaner with nylon bristles in stem. Single use.

K	PC-515	.625"	15.87mm	12"	20/pkg
----------	--------	-------	---------	-----	--------

FLEX ULTRA FINE CLEANER

Polyester cleaner with nylon bristles in stem. Rolls allow you to cut to the size you need, minimizing waste. Single use.

L	PC-606	.236"	6mm	6' roll	1 roll
----------	--------	-------	-----	---------	--------



Instrument Cleaning – Toothbrush-Style

TOOTHBRUSH-STYLE CLEANING BRUSHES

Antimicrobial nylon, brass, or stainless bristles on a sturdy, plastic handle. Designed for external cleaning.

Number	Details	Overall Length	Quantity/pkg
A N-2121	nylon, double-ended, three rows/one row	7"	2
N-2121-50			50
B S-2121	stainless steel, double-ended, three rows/one row	7"	2
S-2121-50			50
C N-2100	nylon, three rows	7"	3
N-2100-50			50
D B-2100	brass, three rows	7"	3
B-2100-50			50
E S-2100	stainless steel, three rows	7"	3
S-2100-50			50
F S-2300	soft stainless steel, three rows	7"	3
S-2300-50			50
G S-2131	soft stainless steel, one row	7"	3
S-2131-50			50
H N-2221	extra rigid nylon, double-ended, three rows/one row	6.89"	2
N-2221-50			50
I N-2200	extra rigid nylon, angled handle	9.25"	2
DOUBLE HEADED TOOTHBRUSH			
J D-2100BN	brass, nylon, three rows	7"	3
D-2100BN-50			50
K D-2100SN	stainless steel, nylon, three rows	7"	3
D-2100SN-50			50







A



B



C



D



E



F



G



H

Instrument Cleaning – General Instrument Cleaning Brushes



GENERAL INSTRUMENT CLEANING BRUSHES

Antimicrobial nylon or stainless steel bristles on a sturdy, plastic handle. Designed for external cleaning.

	Number	Details	Overall Length	Quantity/pkg
A	N-3000	nylon	7"	1
	N-3000-3			3
B	N-4000	nylon	7.25"	1
C	S-4000	stainless steel	7.25"	1

ACETABULAR REAMER BRUSH

Designed for cleaning acetabular bone reamer heads. Rigid, antimicrobial nylon bristles. Double-ended. Round 'pigtail' style brush on one end, straight precision style brush on the other.

	Number	Details	Overall Length	Quantity/pkg
D	N-2000	double-ended	7.725"	1
	N-2000-3			3

GENERAL INSTRUMENT CLEANING BRUSH - STAINLESS STEEL BRISTLES

Stainless steel bristles on a twisted stainless steel handle with loop end.

	Number	Details	Diameter inches	Diameter mm	Bristle Area Length	Length	Quantity/pkg
E	BR-2300	stainless steel	.438"	11.13mm	2"	6"	2

GENERAL INSTRUMENT CLEANING BRUSH - STAINLESS STEEL BRISTLES

Stainless steel bristles on a rigid stainless steel handle. No loop. Designed to reach small, challenging spaces.

	Number	Details	Overall Length	Quantity/pkg
F	BR-04-003	stainless steel	4"	1
	BR-04-003-50			50

DOUBLE-ENDED BRUSH







Soak-proof cleaning brush. Nylon bristles on twisted stainless steel handle with rosebud tip.

	Number	Brush Diameter	Bristle Area Length	Overall Length	Quantity/pkg
G	BR-5450	.350" & 1"	2" & 3.25"	12"	1

FINGER LOOP BRUSH

Nylon bristles on twisted stainless steel handle with rosebud tip. Rigid wire with a finger loop end.

	Number	Brush Diameter	Bristle Area Length	Overall Length	Quantity/pkg
H	BR-5455	1"	3"	10.5"	1



BRUSH TIP DETAILS SHOWN ACTUAL SIZE




Instrument Cleaning – Scope Cleaning Brushes

SCOPE CLEANING BRUSHES

Brush designed for scope cleaning. Nylon bristles attached to flexible, colored plastic tubing. Single-ended brush. Single use.

	Number	Brush Diameter	Bristle Area Length	Length	Color	Quantity/pkg
A	BR-075-18	1.8mm	2cm	75cm	blue	1
	BR-075-18-50					50
B	BR-075-24	2.4mm	2cm	75cm	yellow	1
	BR-075-24-50					50
C	BR-075-32	3.2mm	2cm	75cm	green	1
	BR-075-32-50					50
D	BR-075-42	4.2mm	2cm	75cm	pink	1
	BR-075-42-50					50
E	BR-100-18	1.8mm	2cm	100cm	blue	1
	BR-100-18-50					50
F	BR-100-24	2.4mm	2cm	100cm	yellow	1
	BR-100-24-50					50
G	BR-100-32	3.2mm	2cm	100cm	green	1
	BR-100-32-50					50
H	BR-100-42	4.2mm	2cm	100cm	pink	1
	BR-100-42-50					50
I	BR-150-18	1.8mm	2cm	150cm	blue	1
	BR-150-18-50					50
J	BR-150-24	2.4mm	2cm	150cm	yellow	1
	BR-150-24-50					50
K	BR-150-32	3.2mm	2cm	150cm	green	1
	BR-150-32-50					50
L	BR-150-42	4.2mm	2cm	150cm	pink	1
	BR-150-42-50					50
M	BR-200-18	1.8mm	2cm	200cm	blue	1
	BR-200-18-50					50
N	BR-200-24	2.4mm	2cm	200cm	yellow	1
	BR-200-24-50					50
O	BR-200-32	3.2mm	2cm	200cm	green	1
	BR-200-32-50					50
P	BR-200-42	4.2mm	2cm	200cm	pink	1
	BR-200-42-50					50





BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Scope Cleaning Brushes

SCOPE CLEANING BRUSHES (CONTINUED)

Tapered cleaning brush. Nylon bristles attached to flexible plastic tubing. Single-ended brush. Single use.

	Number	Brush Diameter	Bristle Area Length	Length	Quantity/pkg
A	BR-3308	tapered 3mm to 5mm	1cm	230cm	1
	BR-3308-50				50
B	BR-3309	tapered 3mm to 5mm	1cm	27cm	1
	BR-3309-50				50
C	BR-3310	tapered 3mm to 5mm	2cm	230cm	1
	BR-3310-50				50

DOUBLE-ENDED SCOPE CLEANING BRUSHES

Tapered cleaning brush. Nylon bristles attached to flexible plastic tubing. Double-ended brush. Single use.

D	BR-3320	tapered 3mm to 5mm	1cm	230cm	1
	BR-3320-50				50
E	BR-3330	tapered 3mm to 5mm	2cm	230cm	1
	BR-3330-50				50

DOUBLE-ENDED SCOPE CLEANING BRUSHES

Nylon bristles attached to flexible plastic tubing. 3mm brush on one end, 6mm brush on the other. Single use.


F	BR-3340	3mm & 6mm	1.5cm	230cm	1
	BR-3340-50				50

DELICATE SCOPE CLEANING BRUSHES

Delicate pull-through cleaning brush. Nylon bristles attached to twisted stainless steel wire. Single-ended brush. Single use.

G	BR-3030	1.1mm	1cm	65cm	1
	BR-3030-50				50
	BR-3035	1.1mm	1cm	230cm	1
	BR-3035-50				50
H	BR-3040	2.0mm	1cm	120cm	1
	BR-3040-50				50
I	BR-3410	3.0mm	1cm	120cm	1
	BR-3410-50				50
J	BR-3420	3.0mm	1cm	160cm	1
	BR-3420-50				50





BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Scope Cleaning Brushes



PREMIUM SCOPE CLEANING BRUSHES


Nylon bristles on twisted stainless steel wire or spring coil for hard-to-turn corners. Loop handle at one end and stainless steel beaded tip at bristle end. Reusable.

	Number	Details	Brush Diameter	Bristle Area Length	Length	Quantity/pkg
A	BR-3300	spring coil	tapered 4.0mm to 6.0mm	1cm	246cm	1
B	BR-3305	spring coil	tapered 3.0mm to 5.0mm	1cm	128cm	1
C	BR-3701	twisted wire	2.4mm	6mm	95cm	1
D	BR-3702	spring coil	6.0mm	5mm	95cm	1
E	BR-3703	twisted wire	2.4mm	6mm	165cm	1





SCOPE CLEANING BRUSH

Tapered cleaning brush. Nylon bristles attached to flexible plastic tubing. Single-ended brush. Single use.

	Number	Brush Diameter	Bristle Area Length	Length	Quantity/pkg
F	BR-3711-50 BR-3711-100	tapered 4mm-6mm	1cm	240cm	50 100 (10 bags of 10 brushes)
G	BR-3799-50	Brush Kit, BR-3711 with BR-3201		240cm	50



BRUSH TIP DETAILS SHOWN ACTUAL SIZE



Instrument Cleaning – Scope Cleaning Brushes and Brush Assortments



SUCTION WELL CLEANING BRUSH

Nylon bristles on twisted stainless steel wire with a blue plastic handle. Perfect for cleaning suction wells or other short channels. Single use.

	Number	Brush Diameter	Bristle Area Length	Length	Quantity/pkg
A	BR-3201	11.0mm	4.5cm	3.5"	3
	BR-3201-50				50

DOUBLE-ENDED VALVE CLEANING BRUSH

Nylon bristles on twisted stainless steel handle covered with blue handle in the middle. Two diameters of bristles on one brush for convenient cleaning. Single use.

B	BR-3202	5.0mm & 9.0mm	1cm & 2cm	6"	3
	BR-3202-50				50






FINGER LOOP SCOPE BRUSH

Nylon bristles on twisted stainless steel with convenient finger loop handle. Reusable.

C	BR-3307	10.3mm	4cm	4"	1
---	---------	--------	-----	----	---

SCOPE BRUSH ASSORTMENTS

	Number	Includes one of each of the following brushes:	
D	BR-075-99	BR-075-18, BR-075-24, BR-075-32, BR-075-42	4
E	BR-100-99	BR-100-18, BR-100-24, BR-100-32, BR-100-42	4
F	BR-150-99	BR-150-18, BR-150-24, BR-150-32, BR-150-42	4
G	BR-200-99	BR-200-18, BR-200-24, BR-200-32, BR-200-42	4



BRUSH ASSORTMENTS

Original Brush Assortment

Number Includes one package of each of the following brushes:

- A** BR-4199 BR-18-062, BR-18-087, BR-18-105, BR-18-125,
BR-18-250, FT-18-500, N-2100, N-2121, S-2100

Small Diameter Channel Brush Assortment

Includes one package of each of the following brushes:

- B** BR-4299 BR-08-158, BR-12-079, BR-18-062, BR-18-070, BR-18-087,
BR-18-105, BR-18-125, BR-24-197, BR-24-200

Large Diameter Channel Brush Assortment

Includes one package of each of the following brushes:

- C** BR-4399 BR-06-375, BR-08-630, BR-16-276, BR-16-393, BR-18-250,
BR-18-591, BR-18-880, BR-24-312, BR-24-440

Popular Channel Brush Assortment

Includes one package of each of the following brushes:

- D** BR-4499 BR-06-375, BR-08-158, BR-18-062, BR-18-070, BR-18-087,
BR-18-105, BR-18-125, BR-18-250, BR-24-200

Micro-Ophthalmic Brush Assortment

Includes one package of each of the following brushes:

- E** BR-4599 BR-18-039, BR-18-062, BR-18-070, BR-18-087,
BR-18-105, BR-18-125, BR-18-158, N-2100

Laparoscopic Brush Assortment

Includes one package of each of the following brushes:

- F** BR-4699 BR-12-118, BR-16-197, BR-16-393, BR-3202, N-2121

Endoscopic Brush Assortment

Includes one package of each of the following brushes:

- G** BR-4799 BR-12-079, BR-12-118, BR-12-197, BR-16-197,
BR-16-393, BR-24-197, BR-24-393, BR-3202, N-2121

Scope Brush Assortment

Includes one package of each of the following brushes:

- H** BR-4899 BR-3202, BR-3307, BR-3308, BR-3340

Suction Tube Brush Assortment - French Sizes 3-12

Includes one package of each of the following brushes:

- I** BR-4999 BR-08-092, BR-08-131, BR-08-144, BR-12-079, BR-12-118,
BR-12-158, BR-18-039, BR-18-062, BR-18-070, BR-18-105



A ACTUAL LENGTH IS 31"



B ACTUAL LENGTH IS 36"



C



D



E ACTUAL LENGTH IS 28"



F ACTUAL LENGTH IS 28"

Instrument Cleaning – Sterilizer Cleaning Brushes

STERILIZER CLEANING BRUSH

Brush allows for greater reach and cleaning ability. Bristles are at a forward angle to loosen soil without marring surface.

	Number	Details	Dimensions	Quantity/pkg
A	BR-5000	double-surface brush, polyester bristles	6" x 5" brush head, 31" overall length	1

STERILIZER CLEANING BRUSH HANDLE

White plastic threaded handle for use with BR-5400.

	Number	Details	Overall Length	Quantity/pkg
B	BR-5200	threaded handle, plastic	36"	1

STERILIZER CLEANING BRUSH - SCRUB HEAD

Solid white plastic block eliminates bacteria build-up and won't crack, chip, or sliver. Flare-head design assures thorough cleaning action from any angle. Handle sold separately.

	Number	Details	Dimensions	Quantity/pkg
C	BR-5400	flare-head brush, nylon bristles	5.5" x 9.5" brush head	1

PISTOL-GRIP BRUSH

Soak-proof design with pistol-grip handle. Crimped nylon bristles on a block handle that resists chipping and cracking.

	Number	Bristle Length	Bristle Area	Overall Length	Quantity/pkg
D	BR-4000	1.625"	4" x 4"	8"	1

STERILIZER DRAIN CLEANING BRUSH - WITH PLASTIC HANDLE

White nylon bristles, stainless steel wire and soak-proof plastic handle.

	Number	Details	Dimensions	Quantity/pkg
E	BR-5500	long-handled brush, nylon bristles	1.37" diameter, 28" overall length	1

STERILIZER DRAIN CLEANING BRUSH - WITH STAINLESS STEEL HANDLE

Designed for cleaning out drains, lines, and openings. Stainless steel wire with stiff, crimped nylon bristles. Loop handle.




	Number	Details	Dimensions	Quantity/pkg
F	BR-5600	drain/line brush, nylon bristles	2" diameter, 28" overall length	1

STERILIZER CLEANING BRUSH ASSORTMENT


Assortment includes 1 of each of the following items.

Number	Details
BR-6099	BR-5000, BR-5200, BR-5400, BR-5500, and BR-5600






TWO BARS INCLUDED IN BR-7000.



B




A



SHOWN IN USE. BRUSHES NOT INCLUDED.

C



Instrument Cleaning – Brush Organizer and French Catheter Scale



BRUSH ORGANIZER

Stainless steel organizer with a built-in drip tray. Includes 12 permanent pegs and 2 additional pegged bars (4 pegs per bar). For easy organization and storage display of cleaning brushes.

Number	Details	Quantity
A BR-7000	16" x 35" (1.33' x 2.91') wall-mounted organizer	1/pkg

BRUSH ORGANIZER WITH BRUSH ASSORTMENT

B BR-7100	wall-mounted organizer with Original Brush Assortment
------------------	---

Assortment includes 1 package of the following brushes:

BR-18-062	BR-18-087	BR-18-105
BR-18-125	BR-18-250	FT-18-500
N-2100	N-2121	S-2100

PEGGED BAR FOR BRUSH ORGANIZER

Customize your organizer to fit your needs. Stainless steel bar with 4 pegs. Adhesive-backed for universal placement. Bar measures 1" wide by 5" long.

C BR-7050	pegged bar with adhesive backing	1/pkg
------------------	----------------------------------	-------

FRENCH CATHETER SCALE

Convenient two-sided french catheter scale with built-in 6" ruler.

D MJ-370	plastic card with die-cut holes	4/pkg
-----------------	---------------------------------	-------

